

GE Healthcare



# Compact Airway Modules E-CAiO(V)(X), E-CO(V)(X)

Comprehensive respiratory monitoring for anesthesia and critical care applications

This family of compact airway modules is designed to support respiratory monitoring in anesthesia and critical care areas. Depending on the module version and clinical application needed they provide measurements of airway gases, anesthetic agents with identification, Patient Spirometry and gas exchange/metabolics.

## Features

- Airway gases measured by the sidestream method
- Calculated balance gas values
- All parameter values sampled proximal at the patient's airway with a single gas sampling line, D-lite(+)\* or Pedi-lite(+) flow sensor, along with an additional Spirometry tube
- Anesthesia- and critical care-specific water separation systems to support care area workflow
- Numerical and graphical trends of measured data available together with parameters on the monitor's screen to support multiparameter-based clinical decision making
- Detects end inspiratory and end expiratory occlusions automatically and calculates values for Static Plat, Static PEEPi+e and Static Compliance

## Clinical measurements

- CO<sub>2</sub> and N<sub>2</sub>O – GE infrared technology: Inspired and end-tidal values, CO<sub>2</sub> waveform and respiration rate

- Respiration rate – calculated from the CO<sub>2</sub> waveform
- Anesthetic agents – GE infrared technology
  - Measures and identifies all five agents and two agent mixtures: halothane, enflurane, isoflurane, sevoflurane and desflurane
  - MAC (Minimum Alveolar Concentration)
  - MACage with age, temperature and ambient pressure compensation
- Patient oxygen – GE paramagnetic oxygen (O<sub>2</sub>) technology: Inspired, end-tidal and Fi-Et difference, waveform
- Patient Spirometry – Designed to measure true patient values independent of the ventilator with GE-patented D-lite(+) and Pedi-lite(+) flow sensors and gas samplers at the patient airway
  - Numerical values for airway pressure, minute and tidal volumes, compliance, airway resistance and I:E ratio values, and flow and airway pressure waveforms
  - Continuous measurement of intrinsic, extrinsic and total PEEP
  - Pressure-volume and flow-volume loops
  - Ability to store and print up to six loops
  - Recall saved loops to compare to current loop
  - Module keys to save, print or change loops
- Gas exchange – Direct and continuous measurement
  - Oxygen consumption ( $\dot{V}O_2$ ) and carbon dioxide production ( $\dot{V}CO_2$ )
  - Values for energy expenditure (EE) and respiratory quotient (RQ)



## Technical specifications

### General

Sampling rate 200 ±20 ml/min

Automatic compensation for atmospheric pressure variation (500 to 800 mmHg), temperature and CO<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>O, agent cross effect compensation. Parameter display update interval typically breath-by-breath.

Functional alarms for

- Blocked sample line
- Water trap check
- Water trap replacement
- Low gas sample flow

### Letters in the module name stand for

C = CO<sub>2</sub> and N<sub>2</sub>O

Ai = Anesthetic agents and agent identification

O = Patient O<sub>2</sub>

V = Patient Spirometry

X = Gas exchange

### Non-disturbing gases

Ethanol, acetone, methane, nitrogen, nitric oxide, carbon monoxide, water vapor:

Maximum effect of non-disturbing gases on readings:

CO<sub>2</sub> < 0.2 vol%; N<sub>2</sub>O, O<sub>2</sub> < 2 vol%;

anesthetic agents < 0.15 vol%

### Carbon dioxide (CO<sub>2</sub>)

CO<sub>2</sub> waveform

EtCO <sub>2</sub>	End-tidal CO <sub>2</sub> concentration
FiCO <sub>2</sub>	Inspired CO <sub>2</sub> concentration
Measurement range	0 to 15 vol% (0 to 15 kPa, 0 to 113 mmHg)
Accuracy	±(0.2 vol% + 2% of reading)

GE infrared sensor

Adjustable low and high alarm limits for EtCO<sub>2</sub> or FiCO<sub>2</sub>

### Respiration rate (RR)

Measurement range 4 to 60 breaths/min

Detection criteria 1% variation in CO<sub>2</sub>

Adjustable low and high alarm limits for respiration rate; alarm for apnea

### Patient oxygen (O<sub>2</sub>)

O<sub>2</sub> waveform

FiO <sub>2</sub>	Inspired O <sub>2</sub> concentration
EtO <sub>2</sub>	End-tidal O <sub>2</sub> concentration
FiO <sub>2</sub> -EtO <sub>2</sub>	Inspired-expired difference
Measurement range	0 to 100 vol%
Accuracy	±(1 vol% + 2% of reading)
GE differential paramagnetic sensor	
Adjustable low and high alarm limits for FiO <sub>2</sub> or EtO <sub>2</sub>	

### Nitrous oxide (N<sub>2</sub>O)

FiN <sub>2</sub> O	Inspired N <sub>2</sub> O concentration
EtN <sub>2</sub> O	End-tidal N <sub>2</sub> O concentration
Measurement range	0 to 100 vol%
Accuracy	±(2 vol% + 2% of reading) N <sub>2</sub> O ≤ 85%

Alarm for FiN<sub>2</sub>O >82%

Note: N<sub>2</sub>O is only displayed with CARESCAPE\* ANE and PACU software, and AS/3 and S/5 modular monitors with ANE software.

### Anesthetic agent (AA)

Anesthetic agent waveform

FiAA	Inspired anesthetic agent concentration
EtAA	End-tidal anesthetic agent concentration

MAC or MACage value displayed

Agent mixture detection

Measurement range	
Sevoflurane	0 to 8 vol%
Desflurane	0 to 20 vol%
Isoflurane, enflurane, halothane	0 to 6 vol%
Accuracy	±(0.15 vol% + 5% of reading)

### Agent identification

Identification threshold 0.15 vol%

Adjustable high and low alarm limits for EtAA, FiAA

## Patient Spirometry

---

Pressure-volume loop

Flow-volume loop

Airway pressure and flow waveforms

Adjustable low and high alarm limits for Ppeak, PEEPtot and MVexp

Messages for MVexp << MVinsp and for low volumes

Detection through D-lite or Pedi-lite flow sensor and gas sampler with following specifications:

	<b>D-lite(+)</b>	<b>Pedi-lite(+)</b>
<b>Respiration rate</b>	4 to 35 breaths/min	4 to 50 breaths/min
<b>Tidal volume</b>		
Measurement range	150 to 2000 ml	15 to 300 ml
Accuracy	±6% or 30 ml	±6% or 4 ml
<b>Minute volume</b>		
Measurement range	2 to 20 l/min	0.5 to 5 l/min
<b>Airway pressure</b>		
Measurement range	-20 to +100 cmH <sub>2</sub> O	-20 to +100 cmH <sub>2</sub> O
Accuracy	±1 cmH <sub>2</sub> O	±1 cmH <sub>2</sub> O
Display units	cmH <sub>2</sub> O, mmHg, kPa, mbar, hPa	
<b>Flow</b>		
Measurement range	-100 to +100 l/min	-25 to +25 l/min
<b>I:E</b>		
Measurement range	1:4.5 to 2:1	1:4.5 to 2:1
<b>Compliance</b>		
Measurement range	4 to 100 ml/cmH <sub>2</sub> O	1 to 100 ml/cmH <sub>2</sub> O
<b>Airway resistance</b>		
Measurement range	0 to 40 cmH <sub>2</sub> O/l/s	0 to 40 cmH <sub>2</sub> O/l/s

The presence of xenon or helium in the breathing circuit causes incorrect measurement values.

## Sensor specifications

---

	<b>D-lite(+)</b>	<b>Pedi-lite(+)</b>
Dead space	9.5 ml	2.5 ml
Resistance		
at 30 L/min	0.5 cmH <sub>2</sub> O	
at 10 L/min		1.0 cmH <sub>2</sub> O

## Gas exchange and metabolics

---

$\dot{V}O_2$	Oxygen consumption
$\dot{V}CO_2$	Carbon dioxide production
Measurement range	20 to 999 ml/min
EE	Energy expenditure
Measurement range	0 to 6000 kcal/d 0 to 24000 kJ/d
RQ	Respiratory Quotient ( $\dot{V}CO_2/\dot{V}O_2$ )
Measurement range	0.6 to 1.2 in
Accuracy	Valid for respiration rates 4 to 35 breaths/min
FiO <sub>2</sub> < 65%	± 10% or 10 mL
65% ≤ FiO <sub>2</sub> < 85%	± 15% or 15 mL

Detection through D-lite(+) or Pedi-lite(+) flow sensor and gas sampler. (See the measurement ranges and sensor specifications above.)

Adequate  $\dot{V}O_2/\dot{V}CO_2$  values cannot be measured with leaking airway, with FiO<sub>2</sub> higher than 85%, or when N<sub>2</sub>O or xenon is present in ventilation.

## Monitor compatibility

---

CARESCAPE modular monitors

AS/3, CS/3 and S/5 modular monitors using software S-STD94, S-ARK94, S-ANE97, S-ICU97 or later versions

Displayed data, trends and alarms may vary depending on the host device.

## Environmental specifications

### Operating conditions

---

Temperature	10 to 40°C (50 to 104°F)
Relative humidity	10 to 95% non-condensing

### Storage conditions

---

Temperature	-25 to 70°C (-13 to 158°F)
Relative humidity	10 to 95% non-condensing

## Physical specifications

---

Dimensions (H x W x D)	11.2 x 7.5 x 22.8 cm (4.4 x 3.0 x 9.0 in)
Weight	1.6 kg (3.5 lb)

## Warranty

---

One year

GE Healthcare  
8200 W. Tower Ave.  
Milwaukee, WI 53223  
USA

[www.gehealthcare.com](http://www.gehealthcare.com)

GE Healthcare Finland Oy  
Kuortaneenkatu 2  
00510 Helsinki  
Finland

GE Healthcare  
3/F Building # 1,  
GE Technology Park  
1 Hua Tuo Road  
Shanghai 201203  
China

© 2012 General Electric Company – All rights reserved.  
\* GE, GE Monogram, CARESCAPE and D-lite are  
trademarks of General Electric Company.

General Electric Company reserves the right to make  
changes in specification and features shown herein, or  
discontinue the product described at any time without  
notice or obligation. Contact your GE representative for  
the most current information

Always refer to the user manual that accompanies  
the monitor

GE Healthcare Finland Oy, a General Electric company,  
doing business as GE Healthcare.

GE Healthcare, a division of General Electric Company

