



Allia™ IGS 3, Allia™ IGS 5

Pre-Installation Manual

5871351-8EN
Revision 5
US English

Important...X-Ray Protection



WARNING

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. GE HealthCare will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that GE HealthCare, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

Language Policy

DOC0371395 - Global Language Procedure

PARALAJ-MËRIM (SQ-AL)	<p>Ky manual është i disponueshëm në disa gjuhë.</p> <ul style="list-style-type: none"> Nëse një ofrues shërbimi klientësh kërkon një gjuhë të ndryshme nga ato që mundësohen në Portalin e dokumentacionit të klientit, është përgjegjësia e klientit që të ofrojë shërbime përkthimi. Mos u përpiqni të kryeni shërbime në pajisje, pa lexuar dhe kuptuar paraprakisht manualin e shërbimit. Mosrespektimi i këtij paralajmërimi mund të çojë në lëndim të ofruesit të shërbimit, operatorit ose pacientit si pasojë e goditjes elektrike, mekanike ose një rreziku tjetër.
تحذير (AR-SA)	<p>هذا الدليل متوفر بعدة لغات</p> <ul style="list-style-type: none"> إذا كان مقدم الخدمة التابع للعميل يطلب لغة غير تلك المتوفرة في بوابة توثيق العميل، فإنه يقع على عاتق العميل مسؤولية تقديم خدمات الترجمة لا تحاول صيانة الجهاز ما لم تتم استشارة دليل الخدمة هذا وفهمه قد يؤدي عدم مراعاة هذا التحذير إلى إصابة مقدم الخدمة أو المشغل أو المريض من جراء الصدمات الكهربائية أو المخاطر الميكانيكية أو غيرها من المخاطر
ПРЕДУПРЕЖДЕНИЕ (BG)	<p>Това ръководство е налично на няколко езика.</p> <ul style="list-style-type: none"> Ако доставчикът на услуги на даден клиент изисква език, който е различен от осигурените в портала с документация за клиенти, отговорност на клиента е да предостави преводачески услуги. Не се опитвайте да обслужвате оборудването, освен ако не сте се консултирали с това сервизно ръководство и сте го разбрали. Несъблюдаването на това предупреждение може да доведе до нараняване на предоставящия услугите, оператора или пациента вследствие на токов удар, механична или други опасности.
警告 (ZH-CN)	<p>本手册有多种语言版本。</p> <ul style="list-style-type: none"> 如果客户的服务提供商要求使用 Customer Documentation Portal (客户文档门户) 未提供的其他语言, 则客户有责任提供相应的翻译服务。 请勿尝试检修设备, 除非已明确参考并理解本检修手册。 不遵循此警告可能会导致检修服务提供者、操作员或患者受到触电、机械或其他危害的损伤。
警告 (ZH-HK)	<p>本手冊備有多個語言版本。</p> <ul style="list-style-type: none"> 若客戶的服務提供者所需語言版本不在 Customer Documentation Portal (客戶文件入口網站) 所列語言之中, 客戶需自行負責提供翻譯服務。 除非已查閱並理解本檢修手冊, 否則, 請勿嘗試檢修設備。 不遵循此警告可能會導致服務提供者、操作員或患者因為觸電、機械或其他危險而受傷。
警告 (ZH-TW)	<p>本手冊備有多個語言版本。</p> <ul style="list-style-type: none"> 若客戶的服務提供者所需語言版本不在 Customer Documentation Portal (客戶文件入口網站) 所列語言之中, 客戶需自行負責提供翻譯服務。 除非已查閱並理解本檢修手冊, 否則, 請勿嘗試檢修設備。 不遵循此警告可能會導致服務提供者、操作員或患者因為觸電、機械或其他危險而受傷。
UPOZOR-ENJE (HR)	<p>Ovaj je priručnik dostupan na nekoliko jezika.</p> <ul style="list-style-type: none"> Ako serviser klijenta zahtijeva jezik koji nije jedan od jezika dostupnih na portalu s korisničkom dokumentacijom (Customer Documentation Portal), odgovornost je klijenta pružiti uslugu prevođenja. Nemojte pokušavati servisirati opremu ako niste proučili i razumjeli ovaj servisni priručnik. Nepoštovanje ovog upozorenja može izazvati ozljede servisera, rukovatelja ili pacijenta kao posljedicu strujnog udara, mehaničkih ili drugih opasnosti.

VÝSTRAHA (CS)	<p>Tato příručka je k dispozici v několika jazycích.</p> <ul style="list-style-type: none"> • Pokud zákazníkům poskytovatel služeb vyžaduje jiný jazyk než jazyky, které jsou k dispozici na portálu s uživatelskou dokumentací, je odpovědností zákazníka poskytnout překladatelské služby. • Nepokoušejte se provádět servis zařízení, aniž byste prostudovali tuto servisní příručku a porozuměli jí. • Nedodržení tohoto varování může vést ke zranění poskytovatele služeb, obsluhy nebo pacienta, způsobenému úrazem elektrickým proudem či mechanickým nebo jiným nebezpečím.
ADVARSEL (DA)	<p>Denne vejledning fås på flere sprog.</p> <ul style="list-style-type: none"> • Hvis en kundes tjenesteudbyder kræver et andet sprog end dem, der er til rådighed i Kundedokumentationsportalen, er det kundens ansvar at levere oversættelsestjenester. • Undgå at forsøge at udføre service på udstyret, medmindre du har læst og forstået denne servicevejledning. • Hvis du undlader at overholde denne advarsel, kan det føre til skader på servicemedarbejderen, operatøren eller patienten på grund af elektrisk stød, mekaniske eller andre farer.
WAAR-SCHUWING (NL)	<p>Deze handleiding is in verschillende talen beschikbaar.</p> <ul style="list-style-type: none"> • Als de serviceprovider van een klant een andere taal vereist dan de talen die beschikbaar worden gesteld in het Customer Documentation Portal (Klantdocumentatieportaal), is het de verantwoordelijkheid van de klant om vertaalservices te leveren. • Probeer geen service op de apparatuur uit te voeren zonder de servicehandleiding te hebben gelezen en begrepen. • Het negeren van deze waarschuwing kan leiden tot letsel bij de serviceprovider, de operator of de patiënt door elektrische schokken, mechanische of andere gevaren.
WARNING (EN)	<p>This manual is available in several languages.</p> <ul style="list-style-type: none"> • If a customer's service provider requires a language other than those provided in the Customer Documentation Portal, it is the customer's responsibility to provide translation services. • Do not attempt to service the equipment unless this service manual has been consulted and is understood. • Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.
HOIATUS (ET)	<p>Käesolev juhend on saadaval mitmes keeles.</p> <ul style="list-style-type: none"> • Kui kliendi teenusepakkuja vajab juhendit mõnes muus keeles, mida pole kliendidokumentatsiooni portaalis, on kliendi kohustuseks tõlketeenuste osutamine. • Ärge hakake seda seadet hooldama enne, kui olete käesolevat hooldusjuhendit lugenud ja selle sisu mõistnud. • Selle hoiatuse eiramine võib põhjustada hooldusteenuse pakkuja, operaatorile või patsiendile elektrilöögist, mehhaanilistest või muudest ohtudest tulenevaid vigastusi.
VAROITUS (FI)	<p>Tämä opas on saatavilla useilla kielillä.</p> <ul style="list-style-type: none"> • Jos asiakkaan palveluntarjoaja edellyttää muita kuin asiakkaan asiakirjaportaalissa saatavilla olevia kieliä, käänöspalveluiden tarjoaminen on asiakkaan vastuulla. • Lue huolto-opas huolellisesti ennen laitteen huoltotoimenpiteiden suorittamista. • Tämän varoituksen huomiotta jättäminen voi johtaa huollon suorittajan, laitteen käyttäjän tai potilaan loukkaantumiseen sähköiskun, mekaanisen vaaran tai muun vaaran vuoksi.
ATTENTION (FR)	<p>Ce manuel est disponible en plusieurs langues.</p> <ul style="list-style-type: none"> • Si le prestataire de services d'un client nécessite que le manuel soit rédigé dans une autre langue que celles fournies sur le Portail de Documentation Client, il incombe au client de le faire traduire. • Ne pas essayer d'assurer la maintenance de l'équipement sans avoir au préalable consulté et compris les informations contenues dans ce manuel. • Le non-respect de cet avertissement peut entraîner chez le technicien, l'opérateur ou le patient des blessures dues à des dangers électriques, mécaniques ou autres.

<p>WARNUNG (DE)</p>	<p>Dieses Handbuch ist in mehreren Sprachen erhältlich.</p> <ul style="list-style-type: none"> • Wenn ein Dienstleister des Kunden dieses in einer anderen Sprache als der im Kundendokumentationsportal verfügbaren benötigt, liegt es in der Verantwortung des Kunden, Übersetzungsdienstleistungen zu erbringen. • Wartungsarbeiten am Gerät dürfen nur durchgeführt werden, nachdem dieses Wartungshandbuch gelesen und verstanden wurde. • Andernfalls besteht Verletzungsgefahr für den Dienstleister, Bediener oder Patienten durch Stromschlag, mechanische Gefahren oder andere Gefahren.
<p>ΠΡΟΕΙΔΟΠΟΙ ΗΣΗ (EL)</p>	<p>Αυτό το εγχειρίδιο διατίθεται σε διάφορες γλώσσες.</p> <ul style="list-style-type: none"> • Εάν ο πάροχος υπηρεσιών συντήρησης ενός πελάτη χρειάζεται διαφορετική γλώσσα από αυτές που διατίθενται στο Customer Documentation Portal (Πύλη τεκμηριώσεων πελάτη), ο πελάτης είναι υπεύθυνος για την παροχή υπηρεσιών μετάφρασης. • Μην επιχειρήσετε να εκτελέσετε συντήρηση του εξοπλισμού, εάν δεν έχετε διαβάσει και κατανοήσει το παρόν εγχειρίδιο συντήρησης. • Εάν δεν τηρήσετε αυτήν την προειδοποίηση, μπορεί να προκληθεί τραυματισμός του παρόχου υπηρεσιών συντήρησης, του χειριστή ή του ασθενούς λόγω ηλεκτροπληξίας, μηχανικής βλάβης ή άλλου κινδύνου.
<p>אזהרה (HE)</p>	<p>מדריך זה זמין במספר שפות</p> <ul style="list-style-type: none"> • אם ספק שירות של לקוח זקוק לשפה שאינה מופקת ב Customer Documentation Portal (פורטל) באחריות הלקוח לספק את שירותי התרגום, (תיעוד ללקוחות) • אסור לנסות להעניק שירות לציוד לפני עיון במדריך שירות זה והבנת התוכן שלו • פעולה שלא בהתאם לאזהרה זו עלולה לגרום לפציעה של ספק השירות, המפעיל או המטופל כתוצאה מהתחשמלות, סיכונים מכניים או סיכונים אחרים
<p>FIGYELMEZ- TETÉS (HU)</p>	<p>Ez a kézikönyv több nyelven is rendelkezésre áll.</p> <ul style="list-style-type: none"> • Ha az ügyfél szervizszolgáltatója azoktól eltérő nyelvű kézikönyvet szeretne, mint amelyeket az Ügyféldokumentációs portálon biztosítunk, akkor az ügyfél feladata, hogy gondoskodjon a megfelelő fordításról. • Ne próbálkozzon a berendezés szervizelésével anélkül, hogy a jelen szervizkézikönyvet elolvasta és megértette volna. • Ennek a figyelmeztetésnek a figyelmen kívül hagyása áramütés, mechanikai vagy egyéb veszélyek következtében a szervizszolgáltató, a kezelő vagy a páciens sérülését okozhatja.
<p>ADVÖRUN (IS)</p>	<p>Þessi handbók er ááanleg á mörgum tungumálum.</p> <ul style="list-style-type: none"> • Ef þjónustuaðili viðskiptavinar þarfnast annars tungumáls en þessara tungumála er það á ábyrgð viðskiptavinarins að veita þýðingarþjónustu. • Ekki reyna að þjónusta búnaðinn fyrr en búið er að lesa og skilja þessa þjónustuhandbók. • Sé ekki farið eftir þessari viðvörðun getur það valdið meiðslum á þjónustuaðila, notanda eða sjúklingi af völdum raflosts, vélrænna áverka eða annarar hættu.
<p>PERINGATAN (IN)</p>	<p>Manual ini tersedia dalam beberapa bahasa.</p> <ul style="list-style-type: none"> • Jika penyedia layanan pelanggan membutuhkan bahasa selain dari yang disediakan dalam Portal Dokumentasi Pelanggan, merupakan tanggung jawab pelanggan untuk menyediakan layanan penerjemahan. • Jangan berupaya untuk melakukan servis pada peralatan sebelum menyimpan manual servis dan memahami isinya. • Jika peringatan ini tidak ditaati, ini dapat menyebabkan cedera penyedia layanan, operator, atau pasien, akibat sengatan listrik, bahaya mekanis, atau bahaya lainnya.
<p>AVVERTENZA (IT)</p>	<p>Il presente manuale è disponibile in varie lingue.</p> <ul style="list-style-type: none"> • Qualora un fornitore di servizi del cliente richieda una lingua diversa da quelle fornite nel Portale con la documentazione per il cliente, sarà responsabilità del cliente fornire il servizio di traduzione corrispondente. • Non tentare di riparare l'apparecchiatura se non si è prima consultato e compreso il presente manuale di servizio. • Il mancato rispetto di questa avvertenza può provocare lesioni per il fornitore dei servizi, per l'operatore o per il paziente, a causa di scosse elettriche, meccaniche o altri pericoli.

警告 (JA)	<p>本マニュアルは多言語で提供されています。</p> <ul style="list-style-type: none"> お客様のサービスプロバイダが、お客様ドキュメントポータルページで使用されていない言語を必要とする場合は、お客様の責任で翻訳サービスを提供してください。 機器の保守を行う場合は、必ず本サービスマニュアルを読み理解した上で行ってください。 この警告に従わない場合は、サービスプロバイダー、オペレータ、または患者が、感電、機械的異常、またはその他の有害要因によって負傷する恐れがあります。
경고 (KO)	<p>이 설명서는 여러 언어로 제공됩니다.</p> <ul style="list-style-type: none"> 고객의 서비스 제공자가 고객 문서 포털에 제공된 언어가 아닌 다른 언어를 요구하는 경우 번역 서비스를 제공하는 것은 고객의 책임입니다. 이 서비스 설명서를 참고했고 이해하지 않는 한은 해당 장비를 수리하려고 시도하지 마십시오. 이 경고를 지키지 않으면 감전, 기계상의 위험 또는 다른 위험으로부터 서비스 제공자, 사용자 또는 환자가 다칠 수 있습니다.
BRĪDINĀ- JUMS (LV)	<p>Šī rokasgrāmata ir pieejama vairākās valodās.</p> <ul style="list-style-type: none"> Ja klientu apkalpošanas speciālistam ir nepieciešama cita valoda, kas nav piedāvāta klientu dokumentācijas portālā, klienta pienākums ir nodrošināt tulkošanas pakalpojumus. Nemēģiniet veikt aprikojuma apkopi, kamēr nav izlasīta un izprasta apkopes rokasgrāmata. Ja šis brīdinājums netiek ņemts vērā, pakalpojumu sniedzējs, operators vai pacients var tikt savainots elektriskās strāvas trieciena, mehāniskas vai citas bīstamības rezultātā.
ĮSPĖJIMAS (LT)	<p>Šis vadovas yra išverstas į keletą kalbų.</p> <ul style="list-style-type: none"> Jei kliento paslaugų teikėjui reikalingas vertimas į kitą kalbą, kurios nėra kliento dokumentacijos portale, už vertimo paslaugų suteikimą atsako klientas. Neatlikite įrangos techninės priežiūros, kol neperžiūrėjote ir neišsiaiškinote šio techninės priežiūros vadovo. Nepaisant šio įspėjimo dėl elektros smūgio, mechaninio arba kitokio pavojaus gali būti sužalotas paslaugų teikėjas, operatorius arba pacientas.
TWISSIJA (MT)	<p>Dan il-manwal huwa disponibbli f'diversi lingwi.</p> <ul style="list-style-type: none"> Jekk fornitur tas-servizz ta' klijent ikun jeħtieġ lingwa għajr dawk ipprovduti fil-Portal tad-Dokumentazzjoni tal-Klijent, hija r-responsabbiltà tal-klijent li jipprovidi servizzi ta' traduzzjoni. Tippruvax tagħmel service fuq it-tagħmir sakemm ma jkunx gie kkonsultat u mifhum dan il-manwal għas-service. Jekk wieħed jonqos milli josserva din it-twissija, dan jista' jwassal f'korriment lill-fornitur tas-servizz, lill-operatur jew lill-pazjent minn xokk elettriku, mekkaniku, jew perikli oħra.
ADVARSEL (NO)	<p>Denne håndboken er tilgjengelig på flere språk.</p> <ul style="list-style-type: none"> Hvis en kundes tjenesteleverandør krever et annet språk enn de som finnes i dokumentasjonsportalen for kunder, er det kundens ansvar å levere en oversettelsestjeneste. Ikke prøv å utfør service på utstyret med mindre man har konsultert og forstått servicehåndboken. Om denne advarselen ikke følges kan det føre til skade på tjenesteleverandør, operatør eller pasient fra elektrisk støt, mekanisk eller annen fare.
OSTRZEŻE- NIE (PL)	<p>Niniejszy podręcznik jest dostępny w kilku językach.</p> <ul style="list-style-type: none"> Jeżeli serwisant klienta wymaga języka, który nie został udostępniony w portalu dokumentacji klienta, obowiązkiem klienta jest zapewnienie usług tłumaczeniowych. Nie podejmować prób serwisowania urządzenia bez uprzedniego zapoznania się z niniejszym podręcznikiem serwisowym i zrozumienia jego treści. Nieprzestrzeganie tego ostrzeżenia może spowodować obrażenia u serwisanta, operatora lub pacjenta, spowodowane porażeniem prądem, zagrożeniami mechanicznymi lub innymi.

<p>ATENÇÃO (PT-BR)</p>	<p>Este manual está disponível em vários idiomas.</p> <ul style="list-style-type: none"> • Se o prestador de serviços de um cliente necessitar de um idioma diferente dos fornecidos no Portal da Documentação do Cliente, o fornecimento dos serviços de tradução é de responsabilidade do cliente. • Não tente realizar manutenção do equipamento a menos que o manual de serviço tenha sido consultado e seja entendido. • O não cumprimento deste aviso resultará em lesões ao provedor de serviço, operador ou paciente de choque elétrico, mecânico ou outros riscos.
<p>ATENÇÃO (PT-PT)</p>	<p>Este manual está disponível em vários idiomas.</p> <ul style="list-style-type: none"> • Se o fornecedor de serviços de um cliente necessitar de um idioma diferente dos fornecidos no Portal de Documentação do Cliente, é da responsabilidade do cliente assegurar os serviços de tradução. • Não experimente reparar o equipamento sem primeiro consultar, e compreender, o presente manual de assistência. • O incumprimento deste aviso pode resultar em ferimentos para o técnico de reparação, o operador ou o paciente decorrentes de perigos de eletrocussão, mecânicos ou outros.
<p>ATENȚIE (RO)</p>	<p>Acest manual este disponibil în mai multe limbi.</p> <ul style="list-style-type: none"> • Dacă furnizorul de servicii al unui client necesită o limbă diferită de cele furnizate în Customer Documentation Portal (Portalul cu documentație pentru clienți), este responsabilitatea clientului să furnizeze servicii de traducere. • Nu încercați să efectuați întreținerea echipamentului decât dacă ați consultat și ați înțeles acest manual de service. • Nerespectarea acestei avertizări poate duce la rănirea furnizorului de servicii, a operatorului sau a pacientului din cauza șocurilor electrice, mecanice sau a altor pericole.
<p>ПРЕДУПРЕЖ ДЕНИЕ (RU)</p>	<p>Это руководство доступно на нескольких языках.</p> <ul style="list-style-type: none"> • Если поставщику услуг заказчика требуется языковая версия, отличная от предложенных на портале документации для заказчиков, перевод руководства на необходимый язык осуществляется стороной заказчика. • Не начинайте эксплуатацию оборудования без предварительного надлежащего ознакомления с этим руководством. • Если вы проигнорируете это предупреждение, поставщик услуг, оператор или пациент могут получить механические травмы, травмы вследствие поражения электрическим током или другие увечья.
<p>UPOZOR- ENJE (SR)</p>	<p>Ovaj priručnik je dostupan na nekoliko jezika.</p> <ul style="list-style-type: none"> • Ako korisnikov serviser zahteva neki drugi jezik osim onih koji su dostupni na portalu sa korisničkom dokumentacijom (Customer Documentation Portal), klijent mora da obezbedi prevod. • Nemojte pokušavati da servisirate opremu ako niste proučili i razumeli ovaj priručnik za servisiranje. • Nepoštovanje ovog upozorenja može da izazove povrede servisera, operatera ili pacijenta kao posledicu strujnog udara, mehaničkih ili drugih opasnosti.
<p>UPOZORNE- NIE (SK)</p>	<p>Táto príručka je k dispozícii v niekoľkých jazykoch.</p> <ul style="list-style-type: none"> • Ak poskytovateľ služieb daného zákazníka požaduje jazyk odlišný od jazykov dostupných na portáli s dokumentáciou pre zákazníkov, za prekladateľské služby zodpovedá zákazník. • Nepokúšajte sa vykonávať servis na zariadení, pokiaľ ste si neprečítali a nepochopili pokyny v servisnej príručke. • Nedodržanie tohto varovania môže byť príčinou úrazu poskytovateľa servisu, obsluhy alebo pacienta v dôsledku zásahu elektrickým prúdom alebo v dôsledku mechanických alebo iných nebezpečenstiev.

OPOZORILO (SL)	<p>Ta priročnik je na voljo v več jezikih.</p> <ul style="list-style-type: none"> • Če ponudnik storitev stranke potrebuje priročnik v jeziku, ki ni na voljo na portalu z dokumentacijo stranke, mora stranka zagotoviti prevod. • Opreme ne poskušajte servisirati, če niste prebrali in razumeli tega servisnega priročnika. • V primeru neupoštevanja tega opozorila lahko pride do telesnih poškodb ponudnika storitev, upravljavca ali pacienta zaradi električnega udara, mehanskih ali drugih nevarnosti.
ADVERTENCIA (ES)	<p>Este manual se encuentra disponible en varios idiomas.</p> <ul style="list-style-type: none"> • Si el proveedor de servicios de un cliente requiere un idioma distinto de los proporcionados en el Customer Documentation Portal (Portal de documentación para clientes), es responsabilidad del cliente proporcionar los servicios de traducción. • No intente realizar el mantenimiento del sistema a menos que haya consultado y comprendido este manual de servicio. • El incumplimiento de esta advertencia puede causar lesiones al suministrador de servicios, el operador o el paciente debido a descarga eléctrica, mecánica u otros riesgos.
VARNING (SV)	<p>Denna manual är tillgänglig på flera språk.</p> <ul style="list-style-type: none"> • Om en kunds tjänsteleverantör behöver ett annat språk än de som tillgängliggjorts på portalen för kunddokumentation är det kundens ansvar att erbjuda översättningstjänster. • Försök inte att reparera utrustningen utan att först rådfråga och förstå denna servicehandbok. • Om denna varning inte beaktas kan det leda till skada för tjänsteleverantör, operatör eller patient genom elektrisk stöt, mekaniska eller andra faror.
DİKKAT (TR)	<p>Bu kılavuz birden fazla dile sunulmaktadır.</p> <ul style="list-style-type: none"> • Bir müşterinin servis sağlayıcısı Müşteri Belgeleri Portalı'nda sağlananlardan farklı bir dil talep ederse çeviri hizmeti sağlamak müşterinin sorumluluğundadır. • Bu servis kılavuzuna başvurmadan ve içeriğini anlamadan ekipman üzerinde servis işlemi yapmayı denemeyin. • Bu uyarıya uyulmaması; elektrik çarpması, mekanik tehlikeler veya başka tehlikelerden ötürü servis sağlayıcı, operatör veya hastanın yaralanmasıyla sonuçlanabilir.
ПОПЕРЕДЖЕННЯ (UK)	<p>Цей посібник доступний кількома мовами.</p> <ul style="list-style-type: none"> • Якщо постачальник послуг замовника використовує мову, яку не вказано на порталі з документацією для замовників, послуги з перекладу має забезпечити замовник. • Не починайте роботу з обладнанням без попереднього належного ознайомлення з посібником із використання. • Якщо ви проігноруйте це попередження, постачальник послуг, оператор або пацієнт можуть зазнати механічних травм, ураження електричним струмом або інших тілесних ушкоджень.
CẢNH BÁO (VI)	<p>Tài liệu hướng dẫn này có sẵn ở một số ngôn ngữ.</p> <ul style="list-style-type: none"> • Nếu nhà cung cấp dịch vụ của khách hàng yêu cầu ngôn ngữ khác với ngôn ngữ được cung cấp trong Cổng Thông Tin Tài Liệu Khách Hàng, khách hàng có trách nhiệm cung cấp dịch vụ dịch thuật. • Không cố bảo dưỡng thiết bị trừ khi đã tham khảo và hiểu rõ hướng dẫn sử dụng này. • Việc không chú ý đến cảnh báo này có thể dẫn đến thương tích cho nhà cung cấp dịch vụ, người vận hành hoặc bệnh nhân do điện giật, nguy hiểm cơ học hoặc các mối nguy hiểm khác.

Revision History

Part/Rev	Date	Reason for Change
5871351-8EN Rev A	2021-06	Initial intermediate release of direction 5871351-8EN
5871351-8EN Rev B	2021-09	Second intermediate release of direction 5871351-8EN after review
5871351-8EN Rev 1	2021-11	First release of direction 5871351-8EN Rev 1
5871351-8EN Rev 2	2023-03	Second release of direction 5871351-8EN Rev 2
5871351-8EN Rev 3	2024-02	Third release of direction 5871351-8EN
5871351-8EN Rev 4	2025-01	Fourth release of direction 5871351-8EN (update of Equipment requirements for Gantry height, PDU, LDM and V-Point options, Clinical Access requirement and Power Cables requirements)
5871351-8EN Rev 5	2025-10	Fifth release of direction 5871351-8EN (Electromagnetic Compatibility information update)

Contents

1	General Requirements.....	13
1.1	Objectives & Overview	13
1.1.1	Object and Scope of this manual	13
1.1.2	Pre-Installation Process.....	13
1.2	Customer Responsibilities	14
1.2.1	Responsibilities of the Purchaser/Customer	14
1.2.2	Equipment Classifications	15
1.2.3	Pre-Installation Checklist.....	16
1.3	Delivery Requirements.....	17
1.3.1	Shipping Information.....	17
1.3.1.1	Product Shipping Information	17
1.3.1.2	Detail of Shipping Information.....	20
1.3.2	Tools and Test Equipment	30
1.3.3	Door Size Requirements	33
1.3.4	Route Survey	34
1.4	Product Storage and Handling Requirements	35
1.4.1	Post-delivery warm up period	35
1.4.2	System Storage.....	35
1.4.3	Handling instructions	36
2	Equipment Requirements	39
2.1	System Components	39
2.1.1	Presentation of the 3 Rooms	39
2.1.2	Description of the System.....	40
2.1.2.1	Core System	40
2.1.2.2	Options.....	45
2.1.2.3	Components location and characteristics.....	53
2.1.3	Dimension Drawings	60
2.2	Room Layouts.....	105
2.2.1	Room Dimension Requirements	105
2.2.2	Room Layout Drawings.....	106
2.2.2.1	Exam Room Layout.....	106
2.2.2.2	Control Room Layout.....	113
2.2.2.3	Technical Room Layout	114
2.2.2.4	ECG Device Room Configurations	120
2.2.3	Room Layout Considerations	122
2.2.3.1	Clinical Access	122
2.2.3.2	Peripheral Equipment	122
2.2.3.3	Patient Environment Equipment	123
2.3	Room Structural Requirements	123
2.3.1	General Policy.....	123
2.3.1.1	Baseplates Mounting	123
2.3.1.2	Substructure for Dual Arm suspension Mounting (for Mavig suspension with fixed point dual arm).....	125
2.3.2	Floor Requirements	126
2.3.2.1	Floor requirements when using provided floor anchors	126
2.3.2.2	Floor Construction Requirement	126
2.3.2.3	Hole dimension and preferred location in concrete floor	128

2.3.2.4 Water Pipe Requirements.....	130
2.3.3 Mounting Requirements.....	132
2.3.3.1 Positioner and Table Floor Mounting	132
2.3.3.2 Gantry Table Floor Preparation Kits (GE supplied)	139
2.3.4 Ceiling Requirements.....	139
2.3.4.1 Third Party Monitor suspension (option).....	139
2.3.4.2 Mavig suspension with rails	139
2.3.4.3 MAVIG suspension with fixed point dual arm	141
2.3.5 Wall Requirements.....	143
2.3.5.1 I-Box Installation.....	143
2.3.5.2 (For IGS Control Center) I-Points.....	144
2.3.5.3 Optional Large Display secondary monitor	145
2.3.5.4 V-Point	145
2.4 Seismic.....	146
3 Special Construction Requirements.....	148
3.1 Radiation Protection	148
3.2 Protection against electromagnetic interference hazards.....	148
4 Environmental Requirements	156
4.1 Humidity, Temperature and Altitude	156
4.1.1 Humidity.....	156
4.1.2 Temperature and Altitude	156
4.2 Heat Output.....	157
4.3 Acoustic Specifications	158
4.4 Room Light	158
4.4.1 Requirements for Lighting.....	158
4.4.2 Windows and Curtains	158
4.4.3 Surgical Lights.....	159
5 Electrical Requirements	160
5.1 System Electrical Ratings.....	160
5.2 Power Distribution Schematics	160
5.2.1 System with 1 or 8 kVA UPS	161
5.2.2 System with Fluoro UPS	162
5.2.3 System with Fluoro UPS and IT Electrical Network.....	162
5.2.4 Power distribution for 3rd party monitors on GE 19" monitors suspensions	163
5.3 Cabling Requirements.....	164
5.3.1 Power Source Characteristics.....	164
5.3.1.1 Transformer.....	164
5.3.1.2 Power Cables (Feeders)	165
5.3.2 Fluoro UPS Cables.....	166
5.3.2.1 Fluoro UPS UL Power Cables.....	166
5.3.2.2 Fluoro UPS CE Power Cables.....	166
5.3.2.3 Fluoro UPS Ethernet Cable	167
5.3.3 EPO Cables and additional buttons	167
5.3.3.1 EPO Cable between the MDP and the PDU.....	167
5.3.3.2 Additional EPO buttons	167
5.4 Mains Disconnect Panel	170
5.4.1 General Information.....	170

5.4.1.1	Introduction	170
5.4.1.2	Pre-Installation	170
5.4.1.3	Spare Parts	171
5.4.2	Mandatory Construction Requirements	171
5.4.2.1	Input Power	171
5.4.2.2	Breakers.....	171
5.4.2.3	Terminal Blocks	172
5.4.2.4	Protective Earth	172
5.4.2.5	Indicators	172
5.4.3	Mandatory EPO Requirements	172
5.4.4	Mandatory LOTO Requirements	173
5.4.5	Mandatory Labeling Requirements.....	173
5.4.6	Other Mandatory Requirements.....	173
5.4.7	Preferred Schematics and Components.....	174
5.4.7.1	Recommended CE Schematics	174
5.4.7.2	Minimum Components Specifications for the CE MDP	175
5.4.7.3	Recommended UL Schematics	176
5.4.7.4	Minimum Components Specifications for the UL MDP.....	177
5.4.8	Checklist	178
5.5	External Interfaces.....	179
5.5.1	X-Ray ON lights.....	179
5.5.2	System ON light.....	179
5.5.3	Room lights	180
5.5.4	Room door interlock.....	180
5.5.5	Video distribution	181
5.6	Additional Full UPS.....	182
5.7	Grounding Schematics.....	182
5.8	Third-Party Monitor Suspension Typical Connections.....	185
5.9	System Cable Information	186
5.9.1	Physical Runs	186
5.9.1.1	Physical Run Synoptic	187
5.9.1.2	System Core Matrix	190
5.9.2	MIS (Master Interconnect System).....	200
5.9.3	Cable Channeling.....	201
5.9.3.1	General	201
5.9.3.2	Conduit.....	201
5.9.3.3	Electrical Ducts	201
6	Communication Requirements.....	204
6.1	Network Requirements	204
6.1.1	General Information.....	204
6.1.2	InSite-RSvP Connection Requirements	204
6.1.3	Connection Configuration Parameters	205
6.2	DICOM Requirements	206
6.3	Privacy and Security Configuration.....	206

1 General Requirements

1.1 Objectives & Overview

1.1.1 Object and Scope of this manual

Object and Scope

In this manual:

- The name "Allia™ IGS 5" is used to designate indifferently Allia™ IGS 3 or Allia™ IGS 5. In this case, the requirements are applicable indifferently to both products,
- The term "IGS system" refers to Allia™ IGS 3 and Allia™ IGS 5 product models.
- Allia™ IGS 3 is a product model based on Allia™ IGS 5.

This document is intended as a guide and information resource to properly plan and prepare a site for the installation of an IGS system.

In addition, this document provides references to the pre-installation documents of the various product included with the System.

These documents are intended to assist the Project Manager of Installation (PMI) and the Site Planner in properly preparing a site for the installation of this system.

It provides pre-installation data, such as site preparation prior to the delivery of the System, environmental and electrical requirements and some additional planning aids.

This revision introduces a new Power Safety unit, on among the following would be installed:

- System Interface Cabinet – PDU NPA with 8 kVA UPS could be the legacy UPS 8 kVA Gen1 and the newly UPS 8 kVA Gen2.
- System Interface Cabinet – PDU NPA with 20 kVA Fluoro UPS

The term "UPS 8 kVA" is used to refer to both Gen1 and Gen2 models. The term "FUPS" is used to refer 20 kVA Fluoro UPS.

WARNING



MAKE SURE THE ROOM PREPARATION COMPLIES WITH LOCAL REGULATIONS AS THE PIM IS NOT INTENDED TO REFLECT ALL OF THEM.

Quebec


GE Healthcare is "GE Santé" in Province of Quebec - Canada.

1.1.2 Pre-Installation Process

Complete the checklists in ROOM LAYOUTS, ELECTRICAL REQUIREMENTS, and GENERAL REQUIREMENTS of this manual. They represent an important part of the pre-installation process. The checklists summarize the required preparations and allow to verify the proper completion of the pre-installation procedures.

You will find hereafter a chart of the information flow in the pre-installation process.

Figure 1-1



1.2 Customer Responsibilities

1.2.1 Responsibilities of the Purchaser/Customer

To ensure that the installation of the System meets the purchaser or customer expectations, it is important to determine who will take responsibility for the various items during the system installation process. To help you in determining these responsibilities, review the following checklists with the customer and assign responsibilities as appropriate:

- [Tools and Test Equipment on page 30](#)
- [Pre-Installation Checklist on page 16.](#)

Contract Changes:

Be sure to inform the customer that the cost of any alteration or modification not specified in the sales contract are liable to the customer.

The following equipment must be installed by the Hospital's Contractors, per room drawings:

1. GE-supplied equipment:

- Gantry and Table baseplate with holes drilled (Per supplied template)
- Gantry and Table baseplate grout
- Gantry and Table baseplate floor anchors
- **(For Suspension with rails)** Monitor suspension stationary rails
- **(For Suspension with fixed point Dual Arm)** Substructure for Dual Arm Suspension (S18391MX)



NOTE


Means necessary to anchor of the Substructure for Dual Arm suspension (anchors, bolts, screws, etc.) are not delivered with the kit and should be provided and designed under customer responsibility.


2. Customer supplied equipment:

- MDP (Mains Disconnect Panel)
- Power cables MDP-PDU
- EPO cable MDP-PDU
- Protective Earth cable MDP-PDU
- Power cable MDP-FUPS
- Power cable FUPS-PDU
- **(For Suspension with fixed point Dual Arm)** Means necessary to anchor of the Substructure for Dual Arm suspension (anchors, bolts, screws, etc.)
- Third-Party Monitor suspension

1.2.2 Equipment Classifications

The following equipment classifications are applicable to the product:

Classification category	Equipment classification
Protection against electric shock	<p>Class I</p> <div style="text-align: center;">  <p>WARNING</p> <p>TO AVOID THE RISK OF ELECTRIC SHOCK, THIS EQUIPMENT MUST ONLY BE CONNECTED TO A SUPPLY MAINS WITH PROTECTIVE EARTH.</p> </div>

Classification category	Equipment classification
Degree of protection against electric shock	<p>(For Omega and Innova IQ Table)</p> <p>Type B applied parts</p>  <p>Applied parts complying with the specified requirements of the IEC60601-1 standard to provide protection against electric shock, particularly regarding allowable patient leakage current and patient auxiliary current include:</p> <ul style="list-style-type: none"> • Table Mattress • Table accessories and detachable parts: table head extender, clear-vu arm, head widener with pad/cushion, width extender with pad/cushion, armboard with thick pad/cushion, rail extender and patient restraint strap with cushion. <p>Considered as applied parts: detector cover, removable rails (sleeve).</p>
Degree of protection against harmful ingress of water	<p>Ordinary equipment (enclosed equipment without protection against ingress of water, IPX0) except:</p> <ul style="list-style-type: none"> • Footswitch (protected against the effects of permanent submersion, IPX8). • Innova^{IQ} OR Table, Control Panel, Touch Panel, Table Panning Device (all protected against splashing, IPX4).
Method(s) of sterilization or disinfection recommended by the manufacturer	<ul style="list-style-type: none"> • Sterilization: not applicable • Disinfection: refer to Operator Manual (Chapter Safety and Regulatory, section Disinfection), Recommended disinfecting agents.
Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide	<p>Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide</p> <p>The system does not fulfill the requirements for AP/APG classification (IEC60601-1)</p>
Mode of operation	Continuous operation with intermittent loading
Types of use/reuse	Multiple patient multiple use

NOTICE

THE SYSTEM CAN ONLY BE INSTALLED IN AN ANESTHETIZING LOCATION IF THAT LOCATION IS CLASSIFIED AS OTHER THAN HAZARDOUS AS PER NFPA 70 CLAUSE 517.60.

NOTICE

THE PRODUCT IS NOT CLASSIFIED AS AP, APG (EQUIPMENT NOT SUITABLE FOR USE IN THE PRESENCE OF A FLAMMABLE ANAESTHETIC MIXTURE WITH AIR OR WITH OXYGEN OR NITROUS OXIDE).

1.2.3 Pre-Installation Checklist

Refer to the document *Vascular Site Ready Checklist - DOC2949062* for standard HPM requirements on Room preparation for Vascular Systems installation.

See also the specific preparation requirements for IGS Systems installation given in sections 3, 4 and 5 of the Tab "Installation Prerequisites" in document *IGS System Installation Prerequisites - DOC2024755*.



NOTE

DOC2949062 and DOC2024755 are available from MyWorkshop.

1.3 Delivery Requirements

1.3.1 Shipping Information

1.3.1.1 Product Shipping Information

Refer to the table below. To obtain shipping information for components not specified in the table, refer to the appropriate component Pre-Installation Manual.

Table 1-1

PRODUCT OR COMPONENT	DIMENSIONS			WEIGHT kg (lbs)	METHOD OF SHIPMENT
	Height mm (in)	Width mm (in)	Depth mm (in)		
Gantry (delivered configuration)	2000 (79) (see *)	1230 (48.4)	2820 (111)	1060 (2,340) (see **)	Shipping Dolly. See Figure 1-2 Gantry on Shipping Dolly - Delivered configuration on page 20
	2300 (90.5)	1380 (54.5)	2900 (114)	1200 (2,645) (see **)	Air shipment. See Figure 1-4 Gantry Air Shipment on page 22
Gantry (on-site shortened configuration)	2000 (79)	1160 (46)	2300 (90)	1060 (2,340) (see **)	See Figure 1-3 Gantry on Shipping Dolly - On-site shortened configuration on page 21
C-FRT Cabinet	2200 (87)	1500 (59)	850 (34)	630 (1,388)	On pallet. See Figure 1-5 C-FRT Cabinet Shipment on page 23
NPA PDU	2020 (80)	985 (39)	567 (22)	386 (851)	On pallet. See Figure 1-6 NPA PDU Shipment on page 24
Omega Table Base Assembly	1035 (41)	2060 (81)	840 (33)	585 (1,290)	On pallet. See Figure 1-7 Omega Table Shipment on page 25
Omega Table Top Assembly	220 (9)	3470 (137)	840 (33)	70 (155)	On pallet. See Figure 1-7 Omega Table Shipment on page 25
Innova ^{IQ} OR table Base Assembly and covers	1160 (45.7)	1000 (39.4)	2150 (84.6)	750 (1,653)	On pallet. See Figure 1-8 Innova^{IQ} Table and Innova^{IQ} OR Table Shipment on page 26 and Figure 1-9 Covers Shipment on page 26

Table 1-1 (Table continued)

PRODUCT OR COMPONENT	DIMENSIONS			WEIGHT kg (lbs)	METHOD OF SHIPMENT
	Height mm (in)	Width mm (in)	Depth mm (in)		
X-Ray Tube housing	960 (37.7)	770 (30.3)	710 (28)	113 (250)	On pallet
X-Ray Tube Chiller	1200 (47.2)	555 (21.8)	610 (24)	120 (264.5)	On pallet
Gantry Requisites	-	-	-	-	On pallet
Cables	-	-	-	-	On pallet
Monitor susp. bridge	640 (25.2)	980 (38.6)	3060 (120.5)	210 (445)	On pallet
Monitor susp. rails	380 (15)	300 (12)	5960 (235)	160 (355)	On pallet
Substructure for Dual Arm suspension (4 19" Monitors MAVIG suspension with fixed point dual arm)	330 (13)	1040 (41)	490 (19.3)	70 (154.3)	On pallet, see Figure 1-10 Shipment of Substructure for Dual Arm suspension on page 27
MAVIG suspension with fixed point dual arm for 4 19" Monitors	1860 (73.2)	2150 (84.6)	900 (35.4)	166 (366)	On pallet, see Figure 1-11 MAVIG suspension with fixed point dual arm for 4 19" monitors Shipment on page 28
Large Display monitor (without integrated protective screen)	895 (35.2)	1390 (54.7)	275 (10.8)	45 (99.2)	On pallet
Large Display monitor (with integrated protective screen)	895 (35.2)	1390 (54.7)	275 (10.8)	54.5 (120.2)	On pallet
LD suspension with rails	1100 (43.3)	1100 (43.3)	1850 (72.8)	390 (860)	On pallet
LD suspension with rails 36 m harness	230 (9)	800 (34.5)	800 (34.5)	62 (134)	On pallet
Substructure for Dual Arm suspension (for Mavig suspension with fixed point dual arm for Large Display Monitor)	330 (13)	1040 (41)	490 (19.3)	70 (154.3)	On pallet, see Figure 1-10 Shipment of Substructure for Dual Arm suspension on page 27
Mavig suspension with fixed point dual arm for Large Display Monitor	1860 (73.2)	2150 (84.6)	900 (35.4)	370 (815.7)	On pallet, see Figure 1-12 Large display MAVIG suspension with fixed point dual arm Shipment on page 29
Fluoro UPS UL	2100 (82.7)	890 (35)	1000 (39.4)	561 (1,235)	On pallet
Fluoro UPS CE	1460 (57.5)	540 (21.3)	857 (33.8)	293 (645)	On pallet
8 kVA Gen1 UPS	2100 (82.6)	640 (25.2)	1040 (40.9)	97 (214 lbs)	On pallet
8 kVA Gen2 UPS	700 (27.6)	640 (25.2)	1040 (40.9)	91 (200.6)	On pallet

Table 1-1 (Table continued)

PRODUCT OR COMPONENT	DIMENSIONS			WEIGHT kg (lbs)	METHOD OF SHIPMENT
	Height mm (in)	Width mm (in)	Depth mm (in)		
IGS Control Center not equipped (option for Innova ^{IQ} OR Table)	1700 (66.9)	700 (27.6)	700 (27.6)	46 (101.4)	On pallet
(For IGS Control center) System cables group Cart option	-	-	-	45 (99)	On pallet



NOTE


*: Height can be adjusted: ONLY when necessary on delivery path and IF floor rolling surface is flat and leveled (no obstacle), Dolly can be lowered down by 120-130 mm (it means dolly horizontal bars are at 10 mm from floor surface, to prevent any damage on gantry).

** : If detector is shipped in a separate package, otherwise 29 kg are added.

1.3.1.2 Detail of Shipping Information


1.3.1.2.1 Gantry Shipment Gantry on Shipping Dolly

Figure 1-2 Gantry on Shipping Dolly - Delivered configuration



If the shipping dolly depth is too large you may need to remove the two lift ends and the right and left top handles.

Figure 1-3 Gantry on Shipping Dolly - On-site shortened configuration




Dimensions in mm (in)

Table 1-2 Gantry on Shipping Dolly - On-site shortened configuration

Max width with Gantry mm (in)	Length mm (in)	Height mm (in)
1160 (46)	2300 (90)	2000 (79)

Gantry Air Shipment

Figure 1-4 Gantry Air Shipment



DIMENSIONS IN MM (INCHES)


NOT TO SCALE

1.3.1.2.2 C-FRT Cabinet Shipment

Figure 1-5 C-FRT Cabinet Shipment




Dimensions in mm (in)

 **NOTE** Pallet is delivered as part of C-FRT Cabinet packaging.

1.3.1.2.3 NPA PDU Shipment


Figure 1-6 NPA PDU Shipment



 **NOTE** Pallet is delivered as part of NPA PDU packaging.

1.3.1.2.4 Omega Table Shipment

Figure 1-7 Omega Table Shipment



Dimensions in mm (in)

Item	Description
[1]	Base assembly
[2]	Table top assembly

1.3.1.2.5 Innova^{IQ} Table and Innova^{IQ} OR Table Shipment

Figure 1-8 Innova^{IQ} Table and Innova^{IQ} OR Table Shipment






Figure 1-9 Covers Shipment



1.3.1.2.6 Monitor suspension with fixed point dual arm Substructure for Dual Arm suspension


Figure 1-10 Shipment of Substructure for Dual Arm suspension



Dimensions in mm


MAVIG suspension with fixed point dual arm for 4 19" Monitors

Figure 1-11 MAVIG suspension with fixed point dual arm for 4 19" monitors Shipment



Mavig suspension with fixed point dual arm for Large Display Monitor

Figure 1-12 Large display MAVIG suspension with fixed point dual arm Shipment



Dimensions in mm


1.3.1.2.7 Other Elements Package



NOTE

All OEM parts are shipped inside their original boxes group as needed on pallets.

Figure 1-13 Other Standard Shipping Box



Dimensions in mm (in)

1.3.2 Tools and Test Equipment

Refer to [Table 1-3 on page 30](#). To obtain a list of tools and test equipment for components not specified in [Table 1-3 on page 30](#), refer to the appropriate component Pre-Installation Manual.


Table 1-3


PRODUCT OR COMPONENT	TOOLS	USED FOR	SOURCE	RECEIVED (DATE)
All	Service Engineer's Tool Case	General Use		
LC Positioner	Level, Protractor Type	Positioner Checks		
	Torque Wrench 2 to 20 daN.m (15 ft. lbs. to 150 ft. lbs.)	Positioner Checks		
	1/2 inch Ratchet Wrench	Raise and Lower Positioner shipping dolly		
	Laptop Computer	Positioner Configuration and Calibration		
Omega/Innova ^{IQ} OR Table	Same as for LC Positioner (Service Engineer's Tool Case) Fill in any additional tools or test equipment as required			
	Installation dolly (PN 5265134)	Installation		
Monitor Suspension	Ladders	Installation		
	(For Suspension with rails) XT Lifting Tool (x2) 46-156940G3	Installation		
	(For Suspension with fixed point Dual Arm for 4 19" monitors) • Installation tool and Pelicase (P/N 5758418)	Installation		

Table 1-3 (Table continued)

PRODUCT OR COMPONENT	TOOLS	USED FOR	SOURCE	RECEIVED (DATE)
	<p>(For Suspension with fixed point Dual Arm for LDM)</p> <ul style="list-style-type: none"> • Installation tool and Pelicase (P/N 5758418) • Torque wrench 200 N.m (150 ft.lbs) 	Installation		

Figure 1-14 Omega Table Installation Dolly



NOTE
 Item [A] is delivered mounted on the Omega Table base assembly.

Test Equipments and documentation to store in hospital

All dimensions are in mm and (inch).

- Accessories X-Ray tube cooling unit (to be stored at technical room):
 - 3 boxes : each H x W x D = 210 (8.3) x 400 (15.7) x 310 (12.2)



- Service tools (to be stored in technical room):
 - 2 boxes : each H x W x D = 80 (3.1) x 470 (18.5) x 470 (18.5)



- Suit Cases (to be stored near the exam room as required during calibration procedures):
 - **3D:** H x W x D = 460 (18.1) x 770 (30.3) x 580 (22.8)



- **Positioner:** H x W x D = 90 (3.5) x 340 (13.4) x 260 (10.2)
- **Composite phantom:** H x W x D = 180 (7.1) x 470 (18.5) x 360 (10.2)
- **X-Ray beam accessories:** H x W x D = 130 (5.1) x 450 (17.7) x 380 (15)
- **QAP:** H x W x D = 460 (18.1) x 570 (22.4) x 380 (15)



- Documentation maps (control room or technical room):
 - H x W x D = 330 (13) x 500 (19.7) x 310 (12.2)



1.3.3 Door Size Requirements

Minimum door sizes also apply to hallways and elevators. For additional details, refer to [1.3.1 Shipping Information on page 17](#).

Door Height

The minimum doors height shall be determined to accommodate for the following components:

- The Gantry on its dolly: **2 m** (79 in).
- The C-FRT Cabinet on its pallet: **2.20 m** (87 in).

If the door height is not sufficient, you may need to put the C-FRT Cabinet on its wheels (refer to Dimension Drawings [Figure 2-33 C-FRT Cabinet - Dimensions and CoG on page 80](#) and to *IST0527 - C-FRT Cabinet Installation* in the Service Manual.

WARNING



Adhere to the limit of use described in the Installation Job card.

Door Width

The minimum door width needed (to accommodate the LC shipping dolly) is **1.23 m** (48.4 in) in delivery configuration or **1.16 m** (46 in) in on-site shortened configuration.

NOTE

Door widths are based on a straight-in approach requiring a 2.44 m (96 in) wide corridor. Calculations need to be made for accommodation of equipment through narrower corridors.

Elevator Depth

The minimum elevator depth needed to accommodate the Gantry shipping dolly is:

- **2.82 m** (111 in) with the Gantry shipping dolly in delivery configuration.
- **2.30 m** (90 in) with the Gantry shipping dolly in shortened configuration.


1.3.4 Route Survey

Step One – Sketch

Start preparing Route Survey by sketching a floor plan of the hospital or clinic which will receive the equipment. Include all areas on the delivery route from outside the building to destination. See [Figure 1-15](#) on page 34.

Reference Numbers: Numbers in circles refer to Route Survey data. The Route Survey is a form on which site data are listed (see [Step Two – Survey](#) on page 34).

Figure 1-15



Step Two – Survey

Data concerning the intended delivery route are recorded on the Route Survey in the following pages. Record all loading capacities, corridor widths, door openings, turning radii, flooring materials, elevator sizes, obstructions and so on.

Step Three – Check

Verify equipment can be transported via the route specified in [Step One – Sketch](#) on page 34. Compare Route Survey compiled in [Step Two – Survey](#) on page 34 to equipment specifications in this and other applicable pre-installation directions.

Table 1-6 Storage Requirement

Component	TEMPERATURE		HUMIDITY		PRESSURE	
	MIN	MAX	MIN	MAX	MIN	MAX
All components	+10°C (+50°F)	+40°C (+104°F)	10%	80%	700 hPa	1030 hPa

It is recommended that the temperature for storage does not exceed +25°C (+77°F).

Systems with the Fluoro UPS shall be stored for less than 6 weeks if the storage temperature is above 30°C (86°F), and less than 12 weeks if the storage temperature is above +25°C (+77°F).

Systems with the 1 kVA UPS or 8 kVA UPS shall be stored for less than 14 weeks if the storage temperature is above 30°C (86°F), and less than 25 weeks if the storage temperature is above +25°C (+77°F).

The overall storage time for the system shall be less than 6 months.


Special instructions for the detector:

The detector is very sensitive to temperature and humidity, as irreparable damage to the detector scintillator will occur. If the detector is shipped in a separate package, as defined in [Table 1-6 Storage Requirement on page 36](#), it shall be stored between +10 and +40°C (+50 to +104°F) and less than 80% RH inside its unopened shipping box, the lowest temperature and humidity being preferable. If it is to be stored outside of its shipping box or if the plastic wrapping has been removed, it should be stored at +20°C (68°F) or less and 30% RH or less.

1.4.3 Handling instructions


The packaging of the following components are marked with special handling instructions for transport:

- **Figure 1-16 - Labels on packaging**




- C-FRT Cabinet

Figure 1-17 C-FRT Cabinet - Labels on packaging



- Gantry

Figure 1-18 Gantry - Label on packaging




- Patient Table

Figure 1-19 Patient Table - Label on packaging



- Detector (if shipped separately)

Figure 1-20 Detector - Labels on packaging



2 Equipment Requirements

2.1 System Components

2.1.1 Presentation of the 3 Rooms

The components shall be installed in three different rooms with different constraints: the Exam Room, the Control Room and the Technical Room.

Exam Room

This is where the patient is situated. It contains the table on which the patient is lying, the user interfaces, the gantry, the exam monitors, and accessories.

Control Room

This room contains user interface and control monitors. No intentional or unintentional contact with the patient shall occur with the patient in this area.

Technical Room


This room contains all electronic cabinets. No intentional or unintentional contact with the patient shall occur with the patient in this area. This room shall be separated from the Exam Room and the Control Room, in order to minimize risks of transmission of airborne pathogens. It is recommended to keep the technical room door locked. Its construction should be adapted to minimize ambient noise level; for example the use of glass doors instead of louvered hung doors.

2.1.2 Description of the System

2.1.2.1 Core System

2.1.2.1.1 Gantry

Figure 2-1 Gantry



1. L-arm
2. Pivot
3. X-Ray LED
4. C-Arc
5. Motorized Elevator for the Detector
6. Detector 20 cm, 30 cm or 40 cm
7. X-Ray Tube
8. Collimator
9. X-Ray Tube cover spacer

2.1.2.1.2 Patient Table

Patient Table depends on your System.

Table 2-1 Compatibility matrix for Patient Table

	Allia IGS 320	Allia IGS 330	Allia IGS 520	Allia IGS 530	Allia IGS 540
OMEGA IV Compact Table	Yes	No	Yes	No	No

Table 2-1 Compatibility matrix for Patient Table (Table continued)

	Allia IGS 320	Allia IGS 330	Allia IGS 520	Allia IGS 530	Allia IGS 540
OMEGA V long Table	Yes	Yes	Yes	Yes	Yes
Innova ^{IQ} OR Table	No	No	Yes	Yes	Yes

Figure 2-2 Omega Table (with User Interfaces)



Figure 2-3 Innova^{IQ} OR Table



2.1.2.1.3 User Interfaces

Figure 2-4 User Interfaces



Item	Description
[1]	Control Panel
[2]	Touch Panel
[3]	Table Footswitch
[4]	Table Panning Device (TPD)
[5]	VCIM with DL Keyboard Console
[6]	DLX Keypad
[7]	DL Remote Control

2.1.2.1.4 Accessories for Innova^{IQ} OR Table

Figure 2-5



Item	Description
[1]	Anesthetic Screen Holder
[2]	Round Post clamp
[3]	Sleeve
[4]	Arm Board

Figure 2-6



Item	Description
[1]	Cart
[2]	Head Widener
[3]	Universal Clamp
[4]	Width Extender
[5]	Patient Restraint Strap
[6]	Rail Extender
[7]	User Interface Adaptor Rail

2.1.2.1.5 Monitors

By default:

- Two 19” monitors are provided in the Exam Room:
 - LIVE monitor,
 - REVIEW monitor.

- Two 19” monitors are provided in the Control Room:
 - LIVE monitor,
 - DL monitor.

2.1.2.1.6 Electronic Cabinets and Equipment

The following cabinets and equipments are provided with the system:

- C-FRT Cabinet, which contains the high voltage generator, 2 PCs, IT components and the boards for the Gantry and Table control.
- NPA PDU (Power Distribution Unit) / System Interface Cabinet.
- One UPS among:
 - 1 kVA UPS for system safe shutdown,
 - 8 kVA UPS to maintain all functions except X-Ray acquisitions and Omega table movement during power failures,
 - Fluoro UPS: to complete an exam in fluoroscopy mode during power failures.
- Tube Chiller / Tube Conditioner.
- Detector Conditioner.
- I-Box.
 - **(For System with Omega Table)** mandatory.
 - **(For System with Innova IQ OR Table)** provided only with the IGS Control Center option.

2.1.2.2 Options

2.1.2.2.1 Large Display Monitor (LDM)

The system can integrate a Large Display solution to:

- see images larger at full IQ with greater flexibility in monitor distance in the procedure room.
- display multiple video images simultaneously at different sizes based on stage of workflow.
- conveniently switch operator defined video layouts at different points in procedure workflow.

This option consists in a 55” color monitor and 2 backup 19” monitors in the Exam Room. A second optional Large Display monitor can be provided.

2.1.2.2.2 V-Point

The V-Point option is available only for the systems equipped with the Large Display Option.

The V-Point is a fixed video input for a third party device, located in the Exam Room or in the Control Room. It allows to display the image of this third party device on the LDM. Up to three V-Point can be installed.

The V-Point is compatible with DVI-D (digital only). The maximum supported resolution is 1920 x 1200 60 Hz.

The V-Point can be installed on a wall or on a boom. When installed on the wall, the V-Point **[1]** shall be installed with its box **[2]**.

The V-Point shall not be installed under the table.

Figure 2-7 V-Point Option



The maximum distance between the V-Point and the C-FRT cabinet is 36 m. The diameter of the cable is 20 mm. The routing of the cable shall respect a minimum bending radius of 30 mm when installed on the wall and a minimum dynamic bending of 50 mm radius when mounted on a boom.

2.1.2.2.3 User Interfaces

The user interfaces available on option in the Exam Room are:

- Bolus Handle [1].
- Wireless Footswitch [2].
- **(For System with Omega Table)** Touch Panel 2nd set.

Figure 2-8 Optional User Interfaces



(For System with Innova IQ OR Table) On option, the Control Panel and the Touch Panel can be installed on the IGS Control center. The IGS Control Center is connected to one of the 2 I-Points of the Exam Room.


Figure 2-9 IGS Control Center



The user interfaces available on option in the Control Room are:

- Remote Control Panel connected to the system with the Remote Box.
- In-room AW mouse interface kit (if not in the Exam Room)

Figure 2-10 Remote Box



1



2

Item	Description
[1]	Remote Box for Omega table
[2]	Remote Box for Innova ^{IQ} OR table

2.1.2.2.4 Monitors

According to the subscribed options, up to 9 optional monitors can be installed:

- 1 additional 19" monitor in the Exam Room (Roadmap),
- up to 2 additional 19" monitors in the Control Room (Review and Roadmap),
- up to 6 additional 19" monitors in the Exam Room or in the Control Room.

Table 2-2 Location of 19" monitors (mandatory and optional) - Without LDM

Video Splitter Output	Output 1	Output 2	Output 3	Output 4
Live	Exam Room	Control Room	Exam Room or Control Room	Exam Room or Control Room
Review	Exam Room	Control Room	Exam Room or Control Room	Exam Room or Control Room
Roadmap	Exam Room	Control Room	Exam Room or Control Room	Exam Room or Control Room

Table 2-3 Location of 19" monitors (mandatory and optional) - With LDM

Video Splitter Output	Output 1	Output 2	Output 3	Output 4
Live	Exam Room	Control Room	Exam Room or Control Room	LDM
Review	Exam Room	Control Room	Exam Room or Control Room	LDM
Roadmap	Exam Room	Control Room	Exam Room or Control Room	LDM



NOTE

Text in **bold** for mandatory 19" monitors.

2.1.2.2.5 Monitor Suspensions

GE provides as option several types of suspensions; alternatively, the customer can install the suspension of his choice (third party monitor suspension), provided all requirements in the paragraph [2.1.2.2.5.3 Third Party monitor suspension according to GEHC specifications on page 49](#) are met.

2.1.2.2.5.1 19" monitors suspension (without LDM)

The system can be equipped with a suspension for 4 or 6 19" monitors.

Two types of suspension can be provided:

- suspension with rails: an XT inboard monitor bridge. A monitor frame support receiving 4 or 6 19" monitors (fixed monitor suspension).
- suspension with fixed point Dual Arm: available only for suspension with 4 19" MAVIG monitors.

These suspensions are delivered and installed by GE.

2.1.2.2.5.2 LDM suspension

For the systems with the LDM option, a specific suspension can be provided:

- a suspension with rails
- a suspension with fixed point dual arm.

These suspensions are delivered and installed by GE.

The two backup monitors are mounted on the back of this suspension for faster access in case of failure.

2.1.2.2.5.3 Third Party monitor suspension according to GEHC specifications

The IGS system can be provided with one or several kits to interface third-party suspensions in addition or in replacement of the standard Mavig Suspension usually provided with the system. These kits provide the power and signal connections for the 19" monitors and the Large Display (LD) monitors (8 Megapixels), and for the infra-red receiver and optional dose display.



The association of the IGS system or the purchaser's suspension(s) is not covered by the GEHC product certification.

The overhead monitor suspension shall be installed by strictly following the GEHC installation instructions. GE specifically disclaims any and all liability arising out of or relating to the use or performance of the monitor suspension (including cables), including, without limitation, any liability or claims relating to patient injury, death, or the reliability of such monitors suspension.

The mechanical installation of the third-party suspension and the electrical installation of the third-party monitors are fully under the customer and the installer responsibility. They shall ensure that the third-party suspension and its cables are installed prior to the GE equipment (gantry, table, cabinet) so that the standard GE Service Process can be followed during the system installation. Monitors installation and connections to the GE equipment shall only be made in presence of a GE service representative.

The electrical installation of the third-party monitors is fully under customer and installer responsibility.

The installer is responsible for ensuring that all requirements from this document are met.

It is recommended that the vendor contacts GE Service representative and reviews the site planning details before the suspension is installed. The position of each suspension in the exam room shall be planned following the recommendations from the chapter Room Layout, in order to reduce the risk of collision between any fixed part of the suspension and the gantry or table.

In addition, the location of each suspension in the exam room shall be compatible with the maximum cable length defined in the tables after.



NOTE

GE will not be responsible for any delay in installation if the suspension is not mounted and its cables routed before GE parts arrive on site.

The customer is responsible for providing and replacing any parts of the Third-party Suspension and of the third-party monitors.

Installation requirements:

- The live and roadmap 19" monitors are mandatory in the exam room.
- In order to maintain the IQ performance of the system, only the video cables provided in the kit shall be used to display images on the monitors provided with the system. No extension or additional restpoint is allowed.
- The CAT video cables for the 19" monitors shall respect a minimum bending radius of 24 mm whatever the position of the suspension.
- The CAT video cables for the LDM shall respect a minimum bending radius of 35 mm whatever the position of the suspension.
- Third-party monitors from external sources can be installed on these suspensions but shall not be powered by the system. Only the monitors displaying the images of the system and the AW monitor can be powered by the system.

It is the customer responsibility to ensure that the following requirements are met.

- Each suspension shall not be electrically motorized for up/down motion.

- Each suspension shall comply with the IEC 60601-1 standard and the applicable standards enforced in the country of installation (e.g. when installed in a European Community country, the suspension(s) shall be CE marked). In addition, for North America each suspension shall comply with UL/Canada deviations.
- Each suspension shall be manually adjustable in height and the force to be applied to lift the suspension when fully equipped shall not exceed 200 N in static in the vertical axis, in order to mitigate the risk of patient jammed between table and monitor suspension when the table is lifted up.
- Each suspension shall be installed in order to mitigate the risk of suspension fall on patients and the risk of collisions with the gantry, the table or any other suspension.
- The weight of the monitors and other parts attached to the suspension shall be in accordance with the maximum load supported by the suspension. See [2.1.2.3 Components location and characteristics on page 53](#) for weight and dimensions. For type of fixation:
 - 19” monitors: VESA 100 x 100 mm.
 - Large Display Monitors: VESA 400 x 400 mm.
 - IR Receiver module and Diamentor Display module: VELCRO.
- Each suspension shall be attached to the ceiling in accordance to the manufacturer’s instructions. It shall withstand the maximum suspension load with safety factors in accordance to applicable standards (at least four).
- Each suspension shall be compatible with the chapter Environmental Requirements.
- When the system is installed in an operating room (OR configuration), each suspension shall be compatible with OR environmental constraints.

NOTICE

IT IS MANDATORY TO EXECUTE THE GROUNDING CONTINUITY AND THE LEAKAGE CURRENT MEASUREMENTS AS DEFINED IN CHAPTER 5 AND 6 OF *THIRD-PARTY MONITOR SUSPENSION SERVICE INSTRUCTION FOR INSTALLATION*.

The kits to interface a third-party suspension contain the following cables:

Table 2-4 Kit for the 19" monitors suspension: all cable length 36 m (34.5 m usable length)

From	To	Cable
PDU	Suspension	1 power cable for the monitors powered by the system (the power distribution to each monitor shall be provided by the customer)
PDU	Suspension	1 Protective Earth cables for the suspension (the Protective Earth distribution to each monitor shall be provided by the customer)
C-FRT Cabinet	Monitors	6 RJ45 video cables for the 19” monitors of the system or the AW workstation
C-FRT Cabinet	Infrared receiver	1 cable with D-Sub 9 connector
Dose monitor control device	Optional dose displays	2 cables for 2 optional dose displays
3rd party video source	Monitors	2 VGA cables

Table 2-5 Kit for the LDM suspension: all cable length 36 m (34.5 m usable length)

From	To	Cable
PDU	monitors	3 power cables for the LDM and the 2 backup monitors

Table 2-5 Kit for the LDM suspension: all cable length 36 m (34.5 m usable length) (Table continued)

From	To	Cable
PDU	monitors	3 separate Protective Earth cables for the LDM and the 2 backup monitors
PDU	Suspension	1 Protective Earth cables for the suspension
C-FRT Cabinet	Monitors	4 RJ45 video cables for the LDM 1 RJ 45 video cables for each of the backup monitors
C-FRT Cabinet	Infrared receiver	1 cable with D-Sub 9 connector
Dose monitor control device	Optional dose displays	2 cables for 2 optional dose displays

Table 2-6 Kit for the additional LDM suspension: all cable length 36 m (34.5 m usable length)

From	To	Cable
PDU	monitor	1 power cable
PDU	monitor	1 separate Protective Earth cable
C-FRT Cabinet	Monitor	4 RJ45 video cables

Table 2-7 Kit for 1 additional in-room 19" monitor: all cable length 24 m (22.5 m usable length)

From	To	Cable
PDU	Monitor	1 power cable
PDU	Monitor	1 separate Protective Earth cable
C-FRT Cabinet	Monitor	1 RJ45 video cable

2.1.2.2.6 Advantage Windows workstation (AW)

The AW workstation option is composed of a workstation, 1 or 2 19" flat panel monitors in the Control Room.

Two optional 19" flat panel monitors can be fixed on the Exam Room suspension, or both AW screens can be displayed on the LDM if the option is present.

An optional In-Room AW mouse interface kit can be installed on all systems equipped with workstations versions AW4.6 or higher. It is used to control the workstation through a wireless USB mouse in the Exam Room.

The wireless mouse is not provided in the kit and shall meet the following requirements:

- The shape of the mouse shall be suited for use under a sterile drape.
- It shall have 3 buttons with scroll: the left/right buttons, the up/down scrolling function and the middle scroll button.
- It shall use no specific driver (it shall use generic HID driver).
- The range of the wireless mouse shall be compatible with the distance between the working position and the location of the wireless receiver.

The wireless receiver shall be installed on the wall of the Exam Room, 1.50 m minimum and 5 m maximum from the table, at a suitable height for easy connection of the wireless dongle (1.10 m recommended).

If installed within 1.50 m of the table, the wireless receiver shall be installed inside a box that cannot be opened without a tool. This box shall be provided by the hospital.

The USB port for the wireless receiver can be installed flush on a plaster wall or on a technical sheath, or with a box for concrete walls (the box is provided in the kit):

Plaster wall




Technical sheath



Concrete wall



As an additional option, a Mouse Tray can be attached on the table side rails.



2.1.2.2.7 In-Room third party mouse interface kit installation

An optional In-Room third party mouse interface kit is used to control any third party device located in the control room through a wireless USB mouse in the Exam Room.

Refer to [2.1.2.2.6 Advantage Windows workstation \(AW\) on page 51](#) for the requirements on the wireless mouse and on the installation of the wireless receiver.

2.1.2.2.8 CENTRICITY CA1000

Refer to *Centricity Cardiology CA 1000 V2.0 Preinstallation Guide* in the OEMs of the System Service Manual.

2.1.2.2.9 MacLab

Refer to: Direction 2097992-007, *Mac-Lab/CardioLab/Centricity Cardiology INW Pre-Installation Manual V6.9.6*.

2.1.2.2.10 Injectors

The injectors certified for use with the system are:

Table 2-8

Certified Injectors	Omega Table	Innova ^{IQ} OR Table
MEDRAD Mark 7 (table mount version)	Yes	Yes
MEDRAD Mark 7 (pedestal version)	Yes	Yes
MEDRAD Mark 7 (ceiling mount version)	Yes	Yes
ACIST CVI (pedestal version)	Yes	Yes
ACIST CVI (table mount version)	Yes	No
Medtron Accutron HP model number 836	Yes	Yes
Medtron Accutron HP model number 832	Yes	Yes

Table 2-8 (Table continued)

Certified Injectors	Omega Table	Innova ^{IQ} OR Table
Medtron Accutron HP model number 890	Yes	Yes
Medtron Accutron HP-D model number 833	Yes	Yes



NOTE

For MEDRAD Mark 7 table mount and ceiling mount, rack connected to C-FRT cabinet is located in technical room.

Table accessory rail load considerations:

Each table rail can withstand a load of 40 kg (88 lbs) at 150 mm (5.9") (60 N.m or 44.25 ft/lbs). Therefore, only a light load not exceeding 5 kg (11 lbs) at 100 mm (0.33 ft) can be mounted on the same table rail as the injector: for example IV pole with its accessories, pressure head, and so on. The front table rail is generally used for the User Interfaces, if needed an optional rail can be installed at table foot end of the Omega V table for other options.

The radiation protection and the injector shall never be installed on the same table rail.

2.1.2.2.11 Third-Party on IGS Touch Panel

This option gives the user the ability to remote **Third-Party controls on Allia IGS Touch Panel**.

The **Third-Party on IGS Touch Panel** option requires a **specific Ethernet network topology and cabling** to be configured on site by Service (at system first installation or after).

2.1.2.3 Components location and characteristics

Table 2-9

PRODUCT OR COMPONENT	Number of components in:			NET WEIGHT kg (lbs)	DIMENSIONS WIDTH x DEPTH x HEIGHT mm (in)	LOAD ON THE FLOOR kg/m2 (lb/ft2)
	Exam room	Control room	Technical room			
Gantry	1	-	-	(For System with 20 cm detector) 710 (1,565) (For System with 30 cm detector) 762 (1,680) (For System with 40 cm detector) 786 (1,733)	see Figure 2-11 Gantry - Dimensions and CoG - Side View on page 60 , Figure 2-12 Gantry - Dimensions and CoG - Top view on page 61 and Figure 2-13 Gantry - Dimensions and CoG - Front view on page 62	2678 (549)

Table 2-9 (Table continued)

PRODUCT OR COMPONENT	Number of components in:			NET WEIGHT kg (lbs)	DIMENSIONS WIDTH x DEPTH x HEIGHT mm (in)	LOAD ON THE FLOOR kg/m2 (lb/ft2)
	Exam room	Control room	Technical room			
Omega Table	1	-	-	741.6 (1,635) See NOTE (1)	see Figure 2-14 Omega IV Compact Patient Table - Dimensions and CoG on page 63 and Figure 2-17 Omega V Long Patient Table - Dimensions and CoG on page 66	3065 (628)
Innova ^{IQ} OR Table		--	-	1017 (2,242) See NOTE (1)	see Figure 2-20 Innoval ^{IQ} OR Table - Dimensions on page 69	2260 (463)
Footswitch	1	-	-	Not applicable	Not applicable	Not applicable
Control Panel without clamping	1	-	-	3.5 (6.6)	320 (12.6) x 190 (7.5)	Not applicable
Control Panel with clamping		-	-	4 (8.8)		Not applicable
Touch Panel without clamping	1	-	-	2.5 (4.4)	350 (13.8) x 220 (8.7) x 110 (4.3)	Not applicable
Touch Panel with clamping		-	-	3 (6)		Not applicable
Touch Panel with arm		-	-	4.5 (10)	Not applicable	Not applicable
IGS Control Center (option for Innoval ^{IQ} and Innoval ^{IQ} OR Table)	1	-	-	29 (63.9)	see Figure 2-30 IGS Control Center - Dimensions on page 78	Not applicable
VCIM	-	1	-	0.95 (2.09)	450 (17.7) x 150 (5.9) x 50 (2)	Not applicable
DL Keypad	-	1	-	1.4 (3)	see Figure 2-32 DL Keypad - Dimensions on page 79	Not applicable
DL Monitor	-	1	-	5.6 (12.3)	Not applicable	Not applicable

Table 2-9 (Table continued)

PRODUCT OR COMPONENT	Number of components in:			NET WEIGHT kg (lbs)	DIMENSIONS WIDTH x DEPTH x HEIGHT mm (in)	LOAD ON THE FLOOR kg/m2 (lb/ft2)
	Exam room	Control room	Technical room			
C-FRT Cabinet	-	-	1	556.2 (1,226)	see Figure 2-33 C-FRT Cabinet - Dimensions and CoG on page 80	643 (132)
NPA PDU / System Interface cabinet	-	-	1	291 (642)	see Figure 2-36 NPA PDU Cabinet / System Interface Cabinet - Dimensions and CoG - Side View / Front View / Side View on page 82	847 (173)
Allia Single Plane PDU (NPA PDU) / System Interface cabinet	-	-		290 (639)	see Figure 2-39 Allia Single Plane PDU (NPA PDU) Cabinet / System Interface Cabinet - Dimensions and CoG on page 84	784 (160)
1 kVA UPS	-	-	1	14.8 (32.6)	see Figure 2-42 1 kVA UPS - Dimensions and CoG on page 86	Not applicable
8 kVA UPS Gen1 UPS	-	-		84 (185)	see Figure 2-43 8 kVA Gen1 UPS - Dimensions and CoG on page 87	Not applicable
8 kVA Gen2 UPS	-	-		78.8 (173.7)	see Figure 2-44 Alternate 8 kVA Gen2 UPS - Dimensions and CoG on page 88	Not applicable
Fluoro UPS UL	-	-		530 (1,169)	see Figure 2-45 Fluoro UPS UL - Dimensions and CoG on page 89	975 (200)

Table 2-9 (Table continued)

PRODUCT OR COMPONENT	Number of components in:			NET WEIGHT kg (lbs)	DIMENSIONS WIDTH x DEPTH x HEIGHT mm (in)	LOAD ON THE FLOOR kg/m2 (lb/ft2)
	Exam room	Control room	Technical room			
Fluoro UPS CE	-	-		280 (617)	see Figure 2-46 Fluoro UPS CE - Dimensions and CoG on page 91	1114 (228)
Tube Chiller / Tube Conditioner	-	-	1	120 (265) See NOTE (2)	see Figure 2-47 X-Ray Tube Chiller - Dimensions and CoG on page 92	424 (87)
Detector Conditioner	-	-	1	14.6 (32)	see Figure 2-48 Detector Conditioner - Dimensions and CoG on page 93	Not applicable
(For Omega Table) I-Box	-	-	1	3 (6.6)	see Figure 2-49 I-Box - Dimensions on page 94	Not applicable
OPTIONS						
Monitors						
19" System Monitor without stand	Up to 6	-	-	4.3 (9.5)	405 (15.9) x 61 (2.4) x 334 (13.1)	Not applicable
AW Monitor	Up to 2	-	-	3.1 (6.8)	411 (16.2) x 67 (2.6) x 348 (13.7)	Not applicable
19" System Monitor with stand	-	Up to 3	-	6.1 (13.4)	405 (15.9) x 61 (2.4) x 334 (13.1)	Not applicable
Large Display Monitor (without integrated protective screen)	Up to 2		-	38 (83.8)	1246 (49) x 136 (5.4) x 719 (28.3)	Not applicable
Protective screen option for Large Display Monitor			-	12 (26.5)	1251 (49.3) x 55 (2.2) x 725 (28.6)	Not applicable
Large Display Monitor (with integrated protective screen)			-	47.5 (104.7)	1246 (49) x 136 (5.4) x 719 (28.3)	Not applicable

Table 2-9 (Table continued)

PRODUCT OR COMPONENT	Number of components in:			NET WEIGHT kg (lbs)	DIMENSIONS WIDTH x DEPTH x HEIGHT mm (in)	LOAD ON THE FLOOR kg/m2 (lb/ft2)	
	Exam room	Control room	Technical room				
19" Backup Monitor for Large Display	2	-	-	3.2 (7)	411 (16.2) x 67 (2.6) x 348 (13.7)	Not applicable	
User Interfaces and accessories							
Additional Control Panel	(For Omega Table) Up to 2: 1 in Exam Room and 1 in Control Room (For Innova IQ OR Table) Up to 2: 1 in Exam Room and 1 in Control Room or 2 in Exam Room (Table + IGS Control Center option)			-	Refer to net weight above	Refer to dimensions above	Not applicable
Additional Touch Panel	1	-	-	Refer to net weight above	Refer to dimensions above	Not applicable	
Bolus Handle	1	-	-	Not applicable	Not applicable	Not applicable	
(For InnovaIQ OR Table) I-Box	-	-	1	Refer to net weight above	Refer to dimensions above	Not applicable	
In-room AW mouse interface kit	1		-	Not applicable	Not applicable	Not applicable	
Table Side Cart	1	-	-	Not applicable	Not applicable	Not applicable	
IR Receiver Module	1	-	-	0.3 (0.7)	112 (4.4) x 31 (1.2) x 76 (3)	Not applicable	
Diamentor Display Module	1	-	-	1 (2.2)	210 (8.3) x 146 (5.7) x 58 (2.3)	Not applicable	
ECG Acquisition Device Module							
Physio box	1		-	Not applicable	see Figure 2-50 ECG Acquisition Device Module - Physio Box dimensions (Optional) on page 94	Not applicable	
Suspension							
Precabled 19" LCD monitor suspension with rails for 4 monitors	1	-	-	102 (225)	see Figure 2-51 Suspension with rails for four 19" monitors - Dimensions (Optional) on page 95	Not applicable	

Table 2-9 (Table continued)

PRODUCT OR COMPONENT	Number of components in:			NET WEIGHT kg (lbs)	DIMENSIONS WIDTH x DEPTH x HEIGHT mm (in)	LOAD ON THE FLOOR kg/m2 (lb/ft2)
	Exam room	Control room	Technical room			
Precabled 19" LCD monitor suspension with rails for 6 monitors	1	-	-	115 (254)	see Figure 2-52 Suspension with rails for six 19" monitors - Dimensions (Optional) on page 96	Not applicable
MAVIG suspension with fixed point Dual Arm for 4 19" Monitors	1	-	-	66.9 (147.5)	see Figure 2-57 MAVIG suspension with fixed point dual arm for four 19" monitors - Dimensions (Optional) (1/2) on page 101 and Figure 2-58 MAVIG suspension with fixed point dual arm for four 19" monitors - Dimensions (Optional) (2/2) on page 101	Not applicable
Substructure for Dual Arm suspension (for MAVIG suspension with fixed point Dual Arm for 4 19" Monitors)	1	-	-	58 (128)	see Figure 2-54 Ceiling Plate of Substructure for Dual Arm suspension - Dimensions on page 98	Not applicable
Precabled LD suspension with rails (self weight without monitor and accessories given)	1	-	-	215 (474)	see Figure 2-55 Large Display Suspension with rails - Dimensions (Optional) on page 99	Not applicable
Precabled LD Mavig suspension with fixed point dual arm	1	-	-	190 (419)	see Figure 2-60 MAVIG suspension with fixed point dual arm for Large Display Monitor - Dimensions (Optional) on page 103	Not applicable

Table 2-9 (Table continued)

PRODUCT OR COMPONENT	Number of components in:			NET WEIGHT kg (lbs)	DIMENSIONS WIDTH x DEPTH x HEIGHT mm (in)	LOAD ON THE FLOOR kg/m2 (lb/ft2)
	Exam room	Control room	Technical room			
Substructure for Dual Arm suspension (for LD Mavig suspension with fixed point dual arm)	1	-	-	58 (128)	see Figure 2-54 Ceiling Plate of Substructure for Dual Arm suspension - Dimensions on page 98	Not applicable
Utility						
V-Point	3	-	-	Not applicable	see Figure 2-62 V-Point Box - Dimensions (Optional) on page 105	Not applicable

**NOTE**

(1): Including patient weight. Patient weight considered is:

- 204 kg (450 lbs) for Omega Table,
- 250 kg (551 lbs) for Innova^{IQ} OR Table.

**NOTE**

(2): Dry weight (not filled)

WARNING


THE COMPONENTS IDENTIFIED AS TO BE INSTALLED IN THE TECHNICAL ROOM ARE NOT CERTIFIED FOR USE OUTSIDE OF THIS AREA. IT IS MANDATORY TO INSTALL THEM IN THE TECHNICAL ROOM.

2.1.3 Dimension Drawings

Refer to this section for:

- the dimensional drawings of the system components,
- the location of the components center of gravity,
- Positioner/table relative position drawings.


Figure 2-11 Gantry - Dimensions and CoG - Side View



Dimensions in mm (in)

Item	Description
[1]	Center of Gravity

Figure 2-12 Gantry - Dimensions and CoG - Top view



Dimensions in mm (in)




NOTE

The LC gantry side movement of +95 to -95 is for motorized movement. With manual movement, a range of +100 to -100 is obtained.

Item	Description
[1]	Center of Gravity
*	Recommended
**	Minimum


Figure 2-13 Gantry - Dimensions and CoG - Front view



Dimensions in mm (in)

Item	Description
[1]	Center of Gravity

Figure 2-14 Omega IV Compact Patient Table - Dimensions and CoG




Dimensions in mm (in)

Item	Description
[A]	Side view
[B]	Top view
[1]	Table top at maximum height and patient weight = 204 kg (450 lbs)
[2]	Omega IV compact patient table pivot point
[3]	Center of Gravity with Table top at maximum height and patient weight = 204 kg (450 lbs)
[4]	LCA H1C positioner isocenter
[5]	Table top travel


Item	Description
*	Cardio
**	Neuro

Figure 2-15 Omega IV Compact Patient Table - Plan at Base



Item	Description
[1]	Center of Gravity
[2]	Table Base
[3]	Hole diameter: 29 mm (1.14")


Figure 2-16 Omega IV Compact Patient Table Interference Regions (Allia IGS 520)



Dimensions in mm (in)

Item	Description
[A]	Side view
[B]	Top view
[1]	Omega IV compact patient table pivot point
[2]	Table top travel
*	Cardio
**	Neuro

Figure 2-17 Omega V Long Patient Table - Dimensions and CoG




Dimensions in mm (in)

Item	Description
[A]	Front view (foot end)
[B]	Side view


Item	Description
[C]	Top view
[1]	500 mm (19.5 in) minimum clearance between table end and nearest object
[2]	Table pivot
[3]	Center of Gravity with Table top at maximum height and patient weight = 204 kg (450 lbs)
[4]	Patient coverage = 1980 (77.95) between mechanical stops
[5]	LCA isocenter
[6]	Table top
[7]	Caution: a distance of 700 mm (27.6 in) is mandatory if the head extender is used
*	Wall

Figure 2-18 Omega V Long Patient Table - Plan at Base



Item	Description
[1]	Center of Gravity
[2]	Table Base
[3]	Hole diameter: 29 mm (1.14")

Figure 2-20 Innova^{IQ} OR Table - Dimensions



Dimensions in mm (in)




NOTE

The table dimensions are correct for the position (longi=0, lat=0, lift=0, tilt=0, rot=0)

Item	Description
[A]	Front view (head end)
[B]	Side view
[C]	Top view
[1]	Table pivot

Item	Description
[2]	LCA isocenter

Figure 2-21 Innova^{IQ} OR Table - CoG




Dimensions in mm (in)

Table 2-10

Item	Description
[1]	Center of Gravity

Figure 2-22 Innova^{IQ} OR Table Interference Regions



Dimensions in mm (in)


Item	Description
[1]	Table pivot
[2]	Medical gas box ⁽¹⁾
[3]	Minimum clearance to wall at max rotation
*	Minimum
**	Ceiling or furniture

Item	Description
***	Wall

NOTICE

(1) The medical gas box under the table is not recommended for the surgical configuration.


Figure 2-23 Patient Table side clearance (CPR access)



Dimensions in mm (in)


Item	Description
[1]	Clearance

Figure 2-24 Table Head Extender - Dimensions - Top view / Side view / Front view




Dimensions in mm (in)

Figure 2-25 Gas box outlets Omega IV Table




Item	Description
[A]	Side view
[B]	Plan view
[1]	Positioner isocenter
[2]	See above for distance to isocenter
[3]	Table base
[4]	Table CL
[5]	Table pivot CL
[6]	Alternate area for physio. or med. gases, 12" (305 mm) height restriction in this area
[7]	Recommended area for physio. or med. gases, 21" (533 mm) height restriction in this area
*	Maximum
**	As required

Figure 2-26 Gas box outlets Omega V Table



Item	Description
[A]	Side view
[B]	Plan view
[1]	Positioner isocenter
[2]	See above for distance to isocenter
[3]	Table base
[4]	Table CL
[5]	Table pivot CL
[6]	Recommended area for physio. or med. gases, 12" (305 mm) height restriction in this area
[7]	Recommended area for physio. or med. gases
*	Maximum
**	As required

Figure 2-27 Gas box outlets InnovalQ OR Table




Item	Description
[A]	Side view
[B]	Plan view
[1]	Table base
[2]	Table pivot CL
[3]	Table CL
[4]	Area for physio. or med. gases, 12" (305 mm) height restriction in this area ⁽¹⁾
[5]	Area for physio. or med. gases ⁽¹⁾
*	Maximum

NOTICE


(1) The utility box under the table is not recommended for the surgical configuration.

Figure 2-28 Gantry and Omega IV Compact Patient Table Relative Positions (Allia IGS 520) - Side View




Dimensions in mm (in)

Figure 2-29 Gantry and Omega V Long Patient Table Relative Positions - Side View




Dimensions in mm (in)

Figure 2-30 IGS Control Center - Dimensions




Dimensions in mm (in)

Figure 2-31 19" Desk Mounted Monitor - CoG




Item	Description
[A]	Front Elevation
[B]	Plan at Base
[1]	Center of Gravity

Figure 2-32 DL Keypad - Dimensions



Dimensions in mm (in)


Figure 2-33 C-FRT Cabinet - Dimensions and CoG



Dimensions in mm (in)


Item	Description
[1]	Center of Gravity

Figure 2-34 C-FRT Cabinet - Top Seismic Anchor (Sheet metal S355MC. 1.0976, thickness 5 mm)




Dimensions in mm

Figure 2-35 C-FRT Cabinet - Anchor Bracket (Sheet metal S355MC. 1.0976, thickness 5 mm)




Dimensions in mm

Figure 2-36 NPA PDU Cabinet / System Interface Cabinet - Dimensions and CoG - Side View / Front View / Side View




Dimensions in mm (in)

Figure 2-37 NPA PDU Cabinet - Top bracket (Sheet metal S355MC. 1.0976, thickness 5 mm)




Dimensions in mm

Figure 2-38 NPA PDU Cabinet - Side bracket (Sheet metal S355MC. 1.0976, thickness 5 mm)



Dimensions in mm


Figure 2-39 Allia Single Plane PDU (NPA PDU) Cabinet / System Interface Cabinet - Dimensions and CoG



*: M12x20 Grade 8.8


Dimensions in mm (in)

Figure 2-40 Allia Single Plane PDU (NPA PDU) Cabinet - Top bracket (Sheet metal SPHC, thickness 4 mm)




Dimensions in mm (in)

Figure 2-41 Allia Single Plane PDU (NPA PDU) Cabinet - Side bracket (Sheet metal SPHC, thickness 4 mm)




Dimensions in mm (in)

Figure 2-42 1 kVA UPS - Dimensions and CoG




Dimensions in mm (in)

Figure 2-43 8 kVA Gen1 UPS - Dimensions and CoG




Dimensions in mm (in)

Figure 2-44 Alternate 8 kVA Gen2 UPS - Dimensions and CoG




Dimensions in mm (in)

Figure 2-45 Fluoro UPS UL - Dimensions and CoG




Dimensions in mm (in)




Item	Description
[A]	Front Elevation
[B]	Plan at Base
[1]	Center of Gravity
[2]	Seismic bracket hole diameter: 0.687 inches

Figure 2-46 Fluoro UPS CE - Dimensions and CoG



Dimensions in mm (in)


Figure 2-47 X-Ray Tube Chiller - Dimensions and CoG



Dimensions in mm (in)

Item	Description
[A]	Front Elevation
[B]	Plan at Base
[1]	Center of Gravity


Figure 2-48 Detector Conditioner - Dimensions and CoG



Dimensions in mm (in)


Item	Description
[A]	Front Elevation
[B]	Side Elevation
[1]	Center of Gravity
[2]	Pre-manufactured mounting bracket with 4-3/8" Hilti KB-TZ Expansion Anchors (minimum embedment (h_{ef})=2")

Figure 2-49 I-Box - Dimensions




Dimensions in mm (in)

Figure 2-50 ECG Acquisition Device Module - Physio Box dimensions (Optional)



Dimensions in mm (in)


Figure 2-51 Suspension with rails for four 19" monitors - Dimensions (Optional)



Dimensions in mm (in)

Item	Description
[1]	Optional spacers kit
[2]	XT stationary rail
*	Minimum
**	Maximum


Figure 2-52 Suspension with rails for six 19" monitors - Dimensions (Optional)



Dimensions in mm (in)


Item	Description
[1]	Optional spacers kit
[2]	XT stationary rail
*	Minimum
**	Maximum

Figure 2-53 19" monitors suspension with rails - CoG



Item	Description
[A]	Elevation
[B]	Plan at Ceiling
[1]	Center of Gravity
[2]	Finished Ceiling
[3]	Monitors and Suspension
[4]	Bridge and Dolly
[5]	Ceiling Track (by GE)
[6]	Longitudinal Rail
[7]	4 monitors
*	Maximum
**	Typical

Figure 2-55 Large Display Suspension with rails - Dimensions (Optional)




Dimensions in mm (in)

Item	Description
[1]	Support cable drape rails axis (CPEGE55)
[2]	Optional spacer kit
[3]	XT stationary rail
[4]	HALFEN or UNISTRUT structure

Item	Description
*	Minimum
**	Maximum


Figure 2-56 Large Display Suspension with rails - CoG



Item	Description
[A]	Elevation
[B]	Plan at Ceiling (ceiling mounted)
[1]	Center of Gravity
[2]	Finished Ceiling
[3]	Monitors and Suspension
[4]	Bridge and Dolly
[5]	Longitudinal Rail
[6]	Ceiling Track (by GE)
*	Maximum


Item	Description
**	Typical

Figure 2-57 MAVIG suspension with fixed point dual arm for four 19" monitors - Dimensions (Optional) (1/2)





Dimensions in mm

Figure 2-58 MAVIG suspension with fixed point dual arm for four 19" monitors - Dimensions (Optional) (2/2)



Dimensions in mm

Figure 2-59 MAVIG suspension with fixed point dual arm for four 19" monitors - CoG



Item	Description
[1]	Structural Bracing System
[2]	Structural Plate: 19" x 1" x 1'-7" (A36 MIN)
[3]	Finished Ceiling
[4]	Ceiling Flange Plate: 20 mm THK, S235JR Steel, F _y =52 ksi MIN
[5]	Center of Gravity
[6]	Use 6- MI6 (GR 8.8) threaded rods from adapter plate to support structure
[7]	14.57" diameter bolt circle

Figure 2-60 MAVIG suspension with fixed point dual arm for Large Display Monitor - Dimensions (Optional)






Figure 2-61 MAVIG suspension with fixed point dual arm for Large Display Monitor - CoG



Item	Description
[1]	Structural Support Plate at Support Structure
[2]	Finished Ceiling
[3]	Ceiling Flange Plate: 20 mm THK, S235JR Steel, $F_y=52$ ksi MIN
[4]	Center of Gravity
[5]	Use 6- MI6 (GR 12.9) threaded rods from ceiling flange to support structure
[6]	Structural Plate: 19" x 1" x 1'-7" (A36 MIN)

Figure 2-62 V-Point Box - Dimensions (Optional)



Dimensions in mm (in)

2.2 Room Layouts

2.2.1 Room Dimension Requirements

Table 2-11 Exam Room dimensions

Configuration	Room Length mm (in)		Room Width mm (in)		Ceiling Height mm (in)	
	Recommended	Minimum [A]	Recommended	Minimum	Recommended	Minimum
Omega IV table ⁽¹⁾ [B]	9750 (383.9)	5470 (215.3)	6000 (236.2)	4400 (173.2)	3050 (120.1)	2710 (106.7) with Dual Arm Suspension 2740 (107.9) with rails and Bridge sus- pension
Omega V table ⁽²⁾ [C]		6036 (237.6) (+200 (7.9) from head side with Head Ex- tender [E])				
Innova ^{IQ} OR table ⁽³⁾ [D]		6136 (241.6) (+200 (7.9) from head side with Head Ex- tender [E])				


(1) Applicable for Cardio 1395 mm (54.9 in).

(2) Applicable for both Angio 1278 mm (50.3 in) and Cardio/Neuro 1395 mm (54.9 in).

(3) Applicable for both Angio 1278 mm (50.3 in) and Neuro 1395 mm (54.9 in).

For Head Extender Dimensional drawing, refer to [Dimension Drawings on page 60](#).

Figure 2-63 Minimum Room Length Dimension



2.2.2 Room Layout Drawings

2.2.2.1 Exam Room Layout

NOTICE

Ensure the air outlet positions in the Exam Room is not close to the monitor suspension carriage or rail.

(For System with V-Point) When installed on a wall, the V-Point box should be installed at a suitable height (between 0.80 m and 1.20 m (2.6 ft to 3.9 ft)) from the floor. It should be located near an electrical distribution such as a cable tray or technical sheath, otherwise provide one to route the cables towards the floor or the ceiling. The cable path through the V-Point wall box can be located on one of the four sides of the box or on the back of the box.

The routing of the cable shall respect a minimum bending radius of 30 mm and a minimum dynamic bending of 50 mm radius when mounted on a boom.


Legend of [Figure 2-64 Exam Room Layout - Monitor suspensions rails parallel to patient table \(to side of table\)](#) on page 107, [Figure 2-65 Exam Room Layout - Monitor suspensions rails parallel](#)

to patient table (over table) on page 108, Figure 2-66 Exam Room Layout - Monitor suspensions rails perpendicular to patient table (Gantry side) on page 109, Figure 2-67 Exam Room Layout - Monitor suspensions rails perpendicular to patient table (table side) on page 110 and Figure 2-71 Exam Room Layout - MAVIG suspension with fixed point with dual arm for Large Display Monitor on page 113:

1. Gantry
2. Patient Table
3. Monitor suspension
4. Monitor frame
5. Injector on pedestal

2.2.2.1.1 Exam Room Layout - Monitor Suspension with rails

Figure 2-64 Exam Room Layout - Monitor suspensions rails parallel to patient table (to side of table)



WARNING



Figure 2-64 Exam Room Layout - Monitor suspensions rails parallel to patient table (to side of table) on page 107 is mandatory with hospital having air flow installation.

Figure 2-65 Exam Room Layout - Monitor suspensions rails parallel to patient table (over table)




Figure 2-66 Exam Room Layout - Monitor suspensions rails perpendicular to patient table (Gantry side)





Figure 2-67 Exam Room Layout - Monitor suspensions rails perpendicular to patient table (table side)




2.2.2.1.2 Exam Room Layout - Monitor Suspension with fixed point dual arm

The suspension ceiling fixation shall be determined taking into account at least:

- Clinical need: with an overall radius coverage of 2.21 m for 19" monitor suspension or 2.03 m for Large Display Monitor suspension. Ensure the monitor(s) will be able to reach the position required by medical staff.
- Parking position.
- Ceiling constraints: other component and air flow.
- Cable output and ceiling trap.

Figure 2-68 Exam Room Layout - MAVIG suspension with fixed point with dual arm for four 19" monitors



Item	Description
A	Control Room
B	Exam Room
C	Technical Room
1	Gantry
2	Table
3	4 Monitor Dual Arm Suspension
4	X-ray Rad Shield

The recommended position for the sub-structure ceiling plate installation as in below drawings:

- If the typical procedure of customer is on the patient right, the recommend position refer to [Figure 2-69 on page 112](#).
- If the typical procedure of customer is on the patient left, the recommend position refer to [Figure 2-70 on page 112](#).



NOTE

The longitude edge of the ceiling plate should be parallel to the table top longitude axis.

The distance from center of substructure to the wall should be at least 1300 mm as in the drawing, otherwise, suspension arm could not be rotated full range.

Figure 2-69






Figure 2-70



Item (Figure 2-69 on page 112 and Figure 2-70 on page 112)	Description
A	ISO Center
B	Table Center

Figure 2-71 Exam Room Layout - MAVIG suspension with fixed point with dual arm for Large Display Monitor





2.2.2.2 Control Room Layout

Motion controls installed in remote location from the Table shall be installed at a location where all the positioner axis are visible by the operator.

(For System with InnovalQ OR Table) The remote motion controls shall not be installed on the longitudinal axis of the Table (to avoid any operator visual dead angle due to tilted table top hiding the patient).

(For System with V-Point) The V-Point wall box should be installed on the wall and at a suitable height (between 0.80 m and 1.20 m (2.6 ft to 3.9 ft)) from the floor. It should be near an electrical distribution such as a cable tray or technical sheath, otherwise provide one to route cables towards the floor or the ceiling. Cable path through the V-Point wall box can be located on one of the four sides of the box or on the back of the box. The routing of the cable shall respect a minimum bending radius of 30 mm.

Figure 2-72 Control Room Layout



All measurements are in mm (inches)

2.2.2.3 Technical Room Layout

NOTICE

CONDENSATION MAY OCCUR ON THE OUTLETS AND PIPES OF THE AIR CONDITIONING SYSTEM, THEREFORE, IT IS CRITICAL TO INSTALL THE CABINETS WHERE THERE IS NO RISK OF WATER DROPS FROM THE AIR CONDITIONER.

2.2.2.3.1 General Requirements


It is not allowed to store objects on cabinet top, or to stack cabinets one on another.

In cases 2 cabinets are installed face to face (both sides of the access way), the clearance width shall be at least 1.2 m.

In order to maintain their cooling capacities:

- The Tube Chiller shall be no more than 5 m (16 feet) above or 8 m (26 feet) below the upper position of the X-Ray Tube.

Figure 2-73 Distance between Tube Chiller and X-Ray Tube




NOTE

The highest point of the water network can be 10 m (32 feet) above the floor of the Technical Room where the Coolix-4100 is located (case where the Technical Room is one floor under the Exam Room).

- The Detector Conditioner shall be no more than 3 m in height below the Detector.

Figure 2-74 Technical Room Layout - Configuration 1 kVA UPS



Dimensions in mm (in)

Figure 2-75 Technical Room Layout - Configuration 8 kVA Gen1 UPS






Figure 2-76 Technical Room Layout - Alternate Configuration 8 kVA Gen2 UPS



Dimensions in mm (in)


Figure 2-77 Technical Room Layout - Configuration 20 kVA UPS (Fluoro UPS UL)



Dimensions in mm (in)

Item	Description
[1]	System tools and documentation
[2]	UPS
[3]	C-FRT Cabinet
[4]	PDU
[5]	X-Ray Tube chiller
[6]	Detector chiller
[7]	MDP

Figure 2-78 Technical Room Layout - Configuration 20 kVA UPS (Fluoro UPS CE)



Dimensions in mm (in)

Item	Description
[1]	System tools and documentation
[2]	Fluoro UPS CE
[3]	C-FRT Cabinet
[4]	PDU
[5]	X-Ray Tube chiller
[6]	Detector chiller
[7]	MDP

2.2.2.3.2 Requirements for Equipment Airflow

If the Technical Room is in a dusty environment, it is strongly recommended to install filters on the air inlet of the Technical Room. These filters can cause reduced speed at the air inlet, and the size of the air inlet has therefore to be dimensioned accordingly

The following distances shall be respected to guarantee proper cooling air exhaust.

C-FRT Cabinet:

- The minimum clearance between the ceiling and the top of the C-FRT Cabinet is 30 cm (11.8 in).

NPA PDU:


- The minimum clearance between the ceiling and the top of the NPA PDU is 30 cm (11.8 in).

Tube Chiller:

- The Chiller can operate normally when installed against a wall or another cabinet (no possible air flow) on 1 side. The following clearance must be respected: back and one side: 40 cm minimum, front: 13 cm.

Detector Conditioner:

- The following 50 cm clearance on the sides must be respected.

Figure 2-79 Detector Conditioner – Minimum clearance**Fluoro UPS:**

- Make sure there is a ventilation air flow, preferably ensured by natural air flow, otherwise by enforced ventilation, so that hydrogen concentration is below 1% (according to the Standard IEC 62040-1-2).
- UL UPS: The minimum clearance between the ceiling and the top of the UPS is 40 cm (15.7 in). The left, right or back sides of the UPS can be positioned against the wall or another cabinet.
- CE UPS: The left and right side of the UPS can be positioned against the wall or another cabinet. The following minimum clearance shall be respected:
 - 50 cm (19.7 in) between the top of the UPS and the ceiling
 - 15 cm (5.9 in) between the back of the UPS and the wall
 - 65 cm (25.6 in) from the front of the UPS

2.2.2.3.3 Requirements for Service Access

A free area in front of the following cabinets shall allow to open fully their doors for service access:

- PDU
- C-FRT Cabinet
- MDP
- Fluoro UPS

C-FRT Cabinet:

- A clearance of 80 mm on the lateral sides of the C-FRT Cabinet. It allows the installation of the anchoring brackets and the full opening of the doors for service access

Tube Chiller:

- Minimum 40 cm is recommended for servicing on the left and right side panel. The right side panel is the main side for maintenance. It is recommended to leave this side accessible. The chiller is equipped with wheels that allow to move it during maintenance to allow access on both sides

2.2.2.4 ECG Device Room Configurations

The ECG connection is compatible with an ECG device in the Control Room or in the Exam Room.

The Analog Output Box option is mandatory to provide an analog output connection to the Physio module (If not present, it can be ordered through the following FRUs):

- 2018971-001 16CH ANALOG OUTPUT CPU INTERFACE OPTION


- 2007557-002 KIT ANALOG OUTPUT BOX W/CABLES
- 2010476-001 BOX CARDIOLAB/MACLAB ANALOG OUTPUT

2.2.2.4.1 ECG device in Control Room

Applicable to GE ECG device as MacLab, CardioLab or ComboLab.

In this configuration, the Physio module is installed in the Control Room.


Figure 2-80 ECG device in Control Room - Connection



2.2.2.4.2 ECG device in Exam Room

In this configuration, the Physio module is in the Exam Room.

Figure 2-81 ECG device in Exam Room - Connection



2.2.3 Room Layout Considerations

2.2.3.1 Clinical Access

Make sure that you plan the room with the following clinical access requirements:

- Provide easy access to the patient table. Stretchers and other mobile hospital equipment must reach the table quickly.
- Minimum "trapping zone" safety clearance around the motorized moving parts of the System is 500 mm (19.7 in). Therefore, the minimum distance between the extreme positions of the System, mainly gantry and Omega V/Innova IQ OR tables (or in case there is a table accessory, the extreme position of the accessory) and any static objects (e.g. wall) must be a minimum of 500 mm (19.7 in).

If room size does not allow 500 mm (19.7 in):

- The following text must be added on the room layout proposal reviewed and approved by customer (translated in customer local language):
"Note that your Allia IGS installation in the selected room does not meet the following minimal requirement: 500 mm (19.7 in) required distance between the System (precise which part) and any stationary object.

Therefore the installer must apply a warning label in this area to remind the Operator about entrapment hazard during system motions."
- In any case, table shall not collide with stationary object.
- A standard crushing label - Reference : ISO 7010 Crushing - W019 – must be applied in the area where the 500 mm (19.7 in) clearance is not met.
- Provide sufficient space around the patient table for the unimpeded conduct of CPR (Cardiac Pulmonary Resuscitation). With the table in this position, the table must be capable of rotating $\pm 45^\circ$
- Clinicians at the patient table must be able to communicate with assistants in the control area.
- There must be an unrestricted view of the video monitors and physiological monitoring equipment from the vascular table. Refer to the section Equipment Requirements in the System Pre-Installation Manual.
- Operators in the control area must have easy access to the control console. However, position the controls (including handswitches) so that the operator cannot take exposures while looking around or standing outside the control booth's lead glass window.
- Operators in the control area must have easy access to video recorders, injector programmers, and service and operating manuals.
- Consult customer on the number and location of nonelectrical lines (air, oxygen, vacuum, water, etc.) in the vascular room.
- For systems with the LDM, make sure the backup monitors are easily accessible to view in case of failure of the LDM. For the systems where the backup monitors are mounted at the back of the LDM, plan a clearance so that the monitor can be flipped at 180°.

2.2.3.2 Peripheral Equipment

Consult hospital personnel regarding additional space requirements for the following types of hospital equipment:

- Sinks
- Oxygen stations

- IV apparatus
- Injectors
- Heart monitoring equipment
- Crash cart
- Ultrasound equipment.

2.2.3.3 Patient Environment Equipment

As defined in the IEC60601-1, the patient vicinity is defined as the space within the room 1.83 m (70.7") beyond the perimeter of the table and extending vertically 2.29 m (90.2") above the floor. Only the following components of the system can be installed within the patient vicinity:

- Table and its accessories
- Monitors
- Injector
- Rad-Shield
- User Interfaces
- In-room AW mouse.

2.3 Room Structural Requirements

2.3.1 General Policy

2.3.1.1 Baseplates Mounting

NOTICE

FLOOR, WALLS AND CEILING STRUCTURAL DESIGN THAT MEET MECHANICAL CONSTRAINTS OF SUPPORTING THE SYSTEM, FALL UNDER THE CUSTOMER'S RESPONSIBILITY.

The customer is responsible for the structural analysis and mounting of the base plates. If GE HealthCare is forced to mount the base plate, the Local Customer Team must hire a structural engineer to design and approve the mounting method and provide GE HealthCare with an engineering report.



NOTE

Any flooring material should be removed from the area where the LC Gantry base plate will be mounted. The grout must be directly applied to a properly primed surface, ensure no flooring material exists between the grout pad and the concrete floor. Remove all mastic as possible but do not use a solvent. Solvents prevent good bonding.




NOTE

If the floor is carpeted it must be cut and removed an additional 2 inches (51 mm) must be temporarily lifted out of the way all around if a dam is to be constructed.

NOTICE

The floor slabs on which the equipment is to be installed must have a levelness of 1 mm (0.04 in) per meter (40 in). Position of baseplates and table basement depends on the type of installation. The types of installation are given below.

Figure 2-82 Types of baseplate installation



NOTICE

The Table baseplate or the System baseplate are mandatory to install the Innova^{IQ} OR Table.


The Innova^{IQ} OR Table must never be installed on grade.

The Omega Table can be installed directly on Grade without baseplate. In this case, in floor predrill diameter is 28 mm (1.1 in) and deep is 155 mm (6.1 in).

NOTICE

The gap between the Table Foot bottom end the Gantry Baseplate bottom end shall be lower than 20 mm (0.97 in). Any bigger gap would make the system incompatible with the Vision Applications.

Gap between Table Foot bottom end the Gantry Baseplate bottom



The preferred installation method for the LC Positioner or the Omega tables is through-bolting. The through-bolting method can be used in all seismic zones. If through-bolting cannot be used, use provided floor anchors instead.

2.3.1.2 Substructure for Dual Arm suspension Mounting (for Mavig suspension with fixed point dual arm)

The customer is responsible for the structural analysis and mounting of the Substructure for Dual Arm suspension in the solid ceiling (in case of a MAVIG suspension with fixed point dual arm). If customer requires GEHC to mount the Substructure for Dual Arm suspension, the customer must hire a structural engineer to design and approve the mounting method and provide GEHC with an engineering report.

Figure 2-83 Medium Height Substructure for Dual Arm Suspension and MAVIG Suspension with Fixed Point Dual Arm for four 19" Monitors





Figure 2-84 Medium Height Substructure for Dual Arm Suspension and MAVIG Suspension with Fixed Point Dual Arm for Large Display Monitor



NOTICE

The Substructure for Dual Arm suspension is mandatory to install the MAVIG suspension with fixed point dual arm.

NOTICE

The lower edge of the Substructure for Dual Arm suspension should be the same height as the lower edge of the false ceiling.

2.3.2 Floor Requirements

2.3.2.1 Floor requirements when using provided floor anchors

The maximum pullout force per provided anchor was calculated assuming:


- A concrete compression strength of **17.24 MPa** at 28 days (which is the minimum required compression strength).
- Anchors installed to the required hole depth of **165.1 mm** minimum, and
- Center of anchor hole to concrete edge distance **79.4 mm**.

Make sure to obtain data on compression strength of the concrete before using floor anchors.

2.3.2.2 Floor Construction Requirement

For Slab Type Floor construction, see [Figure 2-85 Through-Bolt Supplied \(Slab Type Floor construction\)](#) on page 126.

Figure 2-85 Through-Bolt Supplied (Slab Type Floor construction)




NOTE

For specific floor preparation procedures, refer to *Floor Pre-Installation Kit Procedures - P/N 5985495*.

For Pan type floor construction, steel channels must be designed by a local structural engineer to span floor joists. See [Figure 2-86 Through-Bolt Supplied \(Pan Type Floor Construction\)](#) on page 127.


Figure 2-86 Through-Bolt Supplied (Pan Type Floor Construction)



NOTE

For detailed information, refer to *Floor Pre-Installation Kit Procedures - P/N 5985495*.

Figure 2-87 Table floor mounting layout



NOTE

For detailed information, refer to *Floor Pre-Installation Kit Procedures - P/N 5985495*.

NOTICE

Prepare the floor such that the table base plate will be flush with the floor finish surface, considering the thickness of the floor finish material.

For alternative table bolts or seismic area, refer to template drawing shown in *Illustration LC Positioner And Table Floor Mounting Template* contained in [2.3.3 Mounting Requirements on page 132](#).

2.3.2.3 Hole dimension and preferred location in concrete floor

In the Exam Room, the LC Positioner is not placed on a computer floor but directly put on concrete floor, the location of the cable access needs to be carefully planned.

Otherwise, if the cable run is located under the concrete floor, the cables will have to come through the floor and in this case you will need two holes, one for the LC Positioner and the other for the patient table.

The diameter of both holes is specified in [Figure 2-88 Hole location in concrete floor on page 129](#).

Figure 2-88 Hole location in concrete floor




Table 2-12 D distance

	ANGIO mm (in)	CARDIO mm (in)	NEURO mm (in)
Omega IV Compact Table	Not applicable	1395 (54.9)	Not applicable
Omega V Long Table	1278 (50.3)	1395 (54.9)	1395 (54.9)

Table 2-12 D distance (Table continued)

	ANGIO mm (in)	CARDIO mm (in)	NEURO mm (in)
Innova ^{IQ} OR Table	1278 (50.3)	1395 (54.9)	1395 (54,9)

NOTICE

Due to the plastic bushing used in the USA to protect cables from the sharp edges of conduits it is necessary to place the cable conduit inside the table cable access opening but the height of the outcoming conduit plus bushing is limited to 1/2 in (12.7 mm).


NOTE

Refer to table *Chemical anchors Pull out efforts and recommendations* in [2.3.3 Mounting Requirements on page 132](#) for pull out effort on each fixation bolts.

2.3.2.4 Water Pipe Requirements

The system uses 4 water pipes for the cooling of the Tube and the Detector. Local regulation may require that electrical cables and water pipes are ran in separate conduits from the Chillers to the Gantry.


Figure 2-89 Water Conduit location with “Above Grade” anchor kits



I.D. = Internal Diameter.
 Suitable External diameter.
 (1) = 9" in US Standard
 225 mm in metric Standard.

Note: Pipe, junction box and duct or conduit are to be supplied and installed by Customer or customer's Contractor.


Figure 2-90 Water Conduit location with “On Grade” anchor kits



NOTE


For flush floor interface, refer to *Floor Pre-Installation Kit Procedures* P/N 5985495.

Figure 2-91 Cable Curvature with “Above Grade” anchor kits



Note: In case of thru-floor cabling, if the electrical cables are coming from the head side, they will need to have a minimum curvature with a minimum radius of 150 mm (6").
In any other cases (i.e. flush floor) no such curvature is allowed.

Figure 2-92 Cable Curvature with “On Grade” anchor kits



NOTE

In case of on grade cabling, because of the minimum curvature constraint of 150 mm (6”), the cable will have to come from the side between the Gantry and the Patient Table.

2.3.3 Mounting Requirements

2.3.3.1 Positioner and Table Floor Mounting

The distances between the Gantry and the Tables are critical for a proper clinical usage. For this reason, GE HealthCare provides two floor mounting templates to ensure these components are properly placed in relation to one another.

Procedures to secure Gantry and Table Floor Mounting Methods are described in *Floor Pre-Installation Kit Procedures* P/N 5985495-199.

Figure 2-93 Gantry Floor Mounting Methods (1/2)




Figure 2-94 Gantry Floor Mounting Methods (2/2)





Figure 2-95 Inner Base Plate (Water Electric Separator) for Above Grade Floor Anchor Kit






Figure 2-96 Cable Conduit For On-Grade Floor Anchor Kit



Dimensions in mm (inches)

ON GRADE




The key hole has been redesigned to permit a 60 mm (2.5 in.) water conduit to separate access for waterpipes used with a MAXIRAY 150.

Push 76 (3.0) hole in center.


Reinsert slug and tack well in one spot hole is for 63.5 (2.5) internal diameter conduit tack well from bottom if possible.

FRONT VIEW



SIDE VIEW

Figure 2-97 Fixing Bolt Overview



NOTE For more details on Floor plate or Base plates, refer to [Figure 2-98 Gantry and table mounting holes](#) on page 138 and illustration *Hole location in concrete floor* in [2.3.2 Floor Requirements](#) on page 126.

- NOTE**
- With any kind of fixation methods (Bolts, Mechanical anchors or Chemical anchors), the mandatory number of holes used is:
 - Gantry base plate: 12 max and 8 min holes used are acceptable
 - Table base plate: 10 max and 8 min holes used are acceptable
 We cannot have 2 consecutive holes omitted.
 - Dedicated to Floor base place, fixation method is only Chemical anchors. The mandatory number of holes used is 22 max and 12 min holes used are acceptable. We can have only 2 consecutive holes omitted.



NOTE

In case of floor plate configuration, the gantry base plate shall be mounted onto the floor plate, but the table base plate shall not be mounted on the floor plate.

Pull out efforts and recommendations about chemical anchors not provided by GE.

The following table provides the recommended chemical anchors for Table/Gantry base plates and for the floor plate ordered locally that they could be used instead of bolts provided by GE.

Table 2-13 Chemical anchors Pull out efforts and recommendations

	Gantry Base Plate	Table Base Plate	Floor Plate (to be ordered locally)	Table Omega
Mark	I on Figure 2-97 Fixing Bolt Overview on page 136	J on Figure 2-97 Fixing Bolt Overview on page 136	K on Figure 2-97 Fixing Bolt Overview on page 136	A on Figure 2-98 Gantry and table mounting holes on page 138
Pull out effort	736 daN per bolt if 12 used and 1992 daN per bolt if 8 used	1120 daN per bolt if 10 used and 2000 daN per bolt if 8 used	272 daN per bolt if 24 used and 2008 daN per bolt if 12 used	4432 daN per bolt with 4 bolts
Number of holes in the plate	12 max (8 min mandatory)	10 max (8 min mandatory)	24 max (12 min mandatory)	4 mandatory
Recommended chemical anchors example 1	Supplier HILTIHVU adhesive capsule + HAS Anchor rod	Supplier HILTIHVU adhesive capsule + HAS Anchor rod	Supplier HILTIHVU adhesive capsule + HAS Anchor rod	Supplier HILTIHVU adhesive capsule + HAS Anchor rod
Threaded rod	M16 A4-70 / 333 131 5/8	M20 A4-70 / 333 135 3/4	M16 A4-70 / 333 131 5/8	M20 A4-70 / 333 135 3/4
Hole diameter in the floor	18 mm (11/16 in)	24 mm (7/8 in)	18 mm (11/16) in	24 mm (7/8 in)
Hole depth in the floor	125 mm (5 in)	170 mm (6-5/8 in)	125 mm (5 in)	170 mm (6-5/8 in)
Minimum floor thickness	180 mm (7 in)	220 mm (8-1/2 in)	180 mm (7 in)	220 mm (8-1/2 in)
Max Tightening Torque	80 N.m (59 ft-lb)	150 N.m (110 ft-lb)	80 N.m (59 ft-lb)	150 N.m (110 ft-lb)



NOTE

The floor plate ordered locally needs to be in steel.

Refer to supplier technical documents for all specification and installation data about chemical anchors.

Figure 2-98 Gantry and table mounting holes




Table 2-14 D distance

	ANGIO mm (in)	CARDIO mm (in)	NEURO mm (in)
Omega IV Compact Table	Not applicable	1395 (54.9)	Not applicable
Omega V Long Table	1278 (50.3)	1395 (54.9)	1395 (54.9)
Innova ^{IQ} OR Table	1278 (50.3)	1395 (54.9)	1395 (54.9)

2.3.3.2 Gantry Table Floor Preparation Kits (GE supplied)


Refer to *Floor Pre-Installation Kit Procedures* - P/N 5985495.

2.3.4 Ceiling Requirements

2.3.4.1 Third Party Monitor suspension (option)

Attention must be paid to the height of suspended elements of the third party monitor suspension, collision must be avoided with the gantry.

Figure 2-99 Gantry and Third Party Monitor Suspension



Dimension in mm

2.3.4.2 Mavig suspension with rails

Aluminum rails support the In-Room Monitor bridge used in the system X-Ray rooms.

Reference:

For additional details on ceiling requirements for stationary rails, refer to Direction 2393190 -1-1EN, *LCD monitor suspension for 2, 3, 4, 6 or 8 monitors* - *Pre-Installation Manual*.

When evaluating ceiling you must take into account the mounting information below.

2.3.4.2.1 Rail Mounting

Attach stationary rails to structural steel with through-bolts in concrete ceilings. Do not use screw anchors in direct tension.

Mount stationary rails directly to the ceiling slab or to flush-mounted unistrut or halfen structure. In higher rooms with false ceiling, mount stationary rails to rigid vertical members hung from ceiling slab.

Securing a supplementary channel to the bottom of the vertical members and mounting the stationary rails to this channel can greatly reduce the number of vertical members.

The stationary rail support structure must be leveled before installation can begin. Do not assume that any support structure is level within specified tolerances, particularly after removing suspensions from an existing room.


2.3.4.2.2 Bolt Specifications (Mavig suspensions)

- The maximum load per bolt will not exceed **1557 N (350 lbs)**.
- Each bolt must not “pull out” or otherwise fail under a vertically downward *dead* load of **6228 N (1400 lbs)**.

2.3.4.2.3 Select Rails (Mavig suspensions)

Monitor suspension rails in two lengths can be selected. Please refer to the Selectable item process or contact the GE representative.

Figure 2-100 Specifications for a typical 19'-0 (5.79 m) inboard stationary rail mounting interface (both rails ceiling mounted), for Mavig suspension



NOTES: 1. NONE CUMULATIVE ERROR.
2. SPACE BETWEEN LAST 2 HOLES MAY BE LESS THAN 26" (66 cm)

Table 2-15 Stationary rail in different length

Rail length cm (in)	Number of holes	A	INBOARD RAILS
472 (186")	8	A1: 7 x 66 cm = 462 cm 7 x 26" = 182"	S18121RC
579 (228")	9	A2: 8 x 66 cm = 528 cm 8 x 26" = 208"	S18121RA

2.3.4.2.4 Cable Support for Monitor Cables

The rails & cable drape are provided with the system, except for the USA, where the Cable Support Kit must be provided locally by the Customer (e.g. CPGE55 from Unistrut).

2.3.4.3 MAVIG suspension with fixed point dual arm

The Substructure for Dual Arm suspension is used to attach the MAVIG suspension with fixed point dual arm to the solid ceiling. It is used as the bridging element between the solid ceiling and the false ceiling for the installation and the fixation of the suspension.

Also, it provides a hooking point required for the installation and the replacement of the Large Display Monitor.

The Substructure for Dual Arm suspension is mandatory to install the MAVIG suspension with fixed point dual arm for Non-seismic Zones. For Seismic Zone installations, refer to Structural Engineer for appropriate design of the structure for installing the MAVIG suspension system.

For standard site configurations, the distance between the ceiling and the lower edge of the false ceiling should be in a range of minimum 175 mm and maximum 610 mm.

NOTICE

If the distance between the ceiling and the lower edge of the false ceiling is more than 610 mm, Long variation of the Substructure for Dual Arm suspension solution could be proposed by MAVIG.

If the distance between the ceiling and the false ceiling is less than 175 mm, then the middle plate is not installed. Refer to [Table 2-16 on page 141](#).

Table 2-16

Distance between ceiling and false ceiling	Configuration of the Substructure for Dual Arm suspension	Item and Description
Minimum is 175 mm and maximum is 610 mm		<p>[1]: Weight in kg [2]: Ceiling Plate [3]: Middle Plate [4]: Maximum is 155 mm</p>
Less than 175 mm		<p>[5]: Maximum is 175 mm / Minimum is 145 mm</p>



NOTE

When distance between ceiling and false ceiling is less than 175 mm, 145 mm minimum is to secure the insertion of the video connection plate delivered with MAVIG Dual Arm suspension.

The Substructure for Dual Arm suspension is delivered with each system. In the GEHC system catalogue (Pre-Installation item), its purchase number is S18391MX (MAVIG Purchase number GD60D022).

2.3.4.3.1 Substructure for Dual Arm suspension mounting

The length of the Substructure for Dual Arm suspension S18391MX can be adapted to any individual situation (distance between solid ceiling and the lower edge of the false ceiling).

Length calculation and adaptation instruction are provided in the MAVIG substructure assembly instructions DBF0100X (where X may be 1 or higher).


The Substructure for Dual Arm suspension must be fastened to the ceiling using six suitable screws.

These screws must be dimensioned according to the conditions of the ceiling and provided by the customer and must be checked by the structural engineer.

The ceiling plate (Figure 2-54 Ceiling Plate of Substructure for Dual Arm suspension - Dimensions on page 98) must be seated flush to the ceiling in order to ensure optimum load distribution.

The lower edge of the Substructure for Dual Arm suspension (Interface plate **[1]**) should be the same as the height as the lower edge of the false ceiling **[2]**.

Figure 2-101 False ceiling alignment versus interface plate



2.3.4.3.2 Bolt Specifications

The statics of the ceiling structure must be designed for safe, load bearing and transfer of the ceiling substructure with suspension arm system and attached devices plus any other existing load applied to the ceiling. In all cases a civil or structural engineer must be consulted.

The screws used to attach the system to the ceiling are not part of the system.

The Substructure shall be fastened to the ceiling with following specifications:

- The maximum axial load per bolt will not exceed 7210 N.
- The maximum Shear load per bolt will not exceed 957 N.
- The maximum pullout force shall be calculated in accordance with local building codes and it is part of structural analysis done by customer.

2.3.4.3.3 False ceiling specifications

The false ceiling should include an opening around the interface plate to allow service engineers to install and replace the suspension and the Large Display Monitor.


The diameter of the opening should be in the range of:

- 489-620 mm (Figure 2-102 on page 143) for a suspension for a Large Display Monitor.
- 577-620 mm (Figure 2-102 on page 143) for a suspension with 4 19" monitors.

A trapdoor in the false ceiling should be provided to allow service access for cables management after mechanical installation of the suspension.

The distance between the substructure and the trapdoor should be less than 50 cm.

Figure 2-102



[1]: Port diameter of the false ceiling: maximum is 620 mm.


2.3.5 Wall Requirements

2.3.5.1 I-Box Installation

(For Omega Table and for InnovalQ OR table with IGS Control Center) The I-Box is securely fastened to the Technical Room wall with four anchors and fixation screws.

The four anchors and fixation screws are not supplied with the I-Box. The fixations shall be sized to support a load of 15 kg (33 lb). Use anchors and screws recommended for the Technical Room wall material.

Figure 2-103 I-Box fixation points



Dimension in mm (in)

2.3.5.2 (For IGS Control Center) I-Points

The 2 I-Points are the connection points of the IGS Control Center. They shall be fixed to the wall of the Exam Room.

Their position is determined by the customer during the layout consideration.

It is recommended to install the I-Points between 800 mm and 1200 mm from finished floor.

The I-Points can be installed on plaster walls or structural walls.

For plaster walls, the minimum required inner distance between the structural wall and the plaster wall is 70 mm.

For structural walls, the I-Point shall be installed in an additional device (e.g. box) to allow its installation on the wall.

Conduits on the wall are also necessary to route the cables. These devices are not provided by GE and must be designed, calculated and supplied locally.

Means of fixation of the I-Points are not delivered with the system (to be provided by the hospital). I-Points must be securely fastened to the Exam Room walls.


The recommended opening dimensions in the plaster wall or box are:

- Hole diameter 90 mm or hole 100 mm x 70 mm for the I-Points.

I-Point dimensions:

- Length: 140 mm
- Width: 100 mm
- Depth: 75 mm

Figure 2-104




2.3.5.3 Optional Large Display secondary monitor

A hooking point shall be provided in order to lift the monitor on a third-party suspension during installation:

- Hooking point characteristic: It must withstand up to 440 lbs (200 kg).
- Recommended hooking point dimensions:

Figure 2-105 Hooking point dimensions




2.3.5.4 V-Point

The V-Point wall box is attached on the wall with four screws and four flat washers. The V-Point on a boom is attached with two screws.



NOTE

The V-Point screws and washers are not provided with the kit. They should be provided under customer responsibility.

Figure 2-106 V-Point Box attached on the wall

2.4 Seismic

Seismic areas

Consider local seismic requirements when planning cabinet mounting.

Consult seismic expert to determine which mounting method is appropriate for the seismic region. Seismic requirements are determined and specified by the hospital/ Design Professional of record and may require approval by the specific state or country agency. Additional reinforcement in the walls may be required by specific seismic areas.

Contact your local GE Installation Program Manager to obtain the latest seismic calculations per the California Building Code (CBC) and the International Building Code (IBC).

The C-FRT Cabinet and the NPA PDU must be securely fastened to the wall and with their seismic kit to prevent them from tipping.

The C-FRT Cabinet, the Detector Conditioner and the Fluoro UPS are each provided with their own seismic kits, excluding the bolts, that shall be provided locally by the customer.

The following seismic kits can be ordered separately (on option):

- NPA PDU seismic kit: S18761PR
- Monitor Flat Panel seismic kit: 5561139
- VICM seismic kit: 2365510
- 1 or 8 kVA UPS seismic kit: E4502YB.

The seismic kit for the Tube Chiller is provided locally by the customer.

(For LDM Suspension with fixed point Dual Arm) :

CAUTION



The standard substructure (MAVIG GD60D022) should not be used with system in seismic zone.

Contact MAVIG or Local contractor to design and supply specific substructure including M12 threaded holes requirement (see below).

Four M12 threaded holes with hooking point are required for the installation of the dual arm suspension, the installation and replacement of the Large Display Monitor. The structural support plate ([Figure 2-61 MAVIG suspension with fixed point dual arm for Large Display Monitor - CoG on page 104](#)) should include these 4 x M12 threaded holes.

For the threaded holes positioning on the structural support plate refer to [Figure 2-54 Ceiling Plate of Substructure for Dual Arm suspension - Dimensions on page 98](#).

Center of Gravity

Refer to [2.1.3 Dimension Drawings on page 60](#) for location of Center of Gravity (CoG) of the system components.

3 Special Construction Requirements

3.1 Radiation Protection

Because X-ray equipment produces radiation, special precautions may be needed or special site modifications may be required. GEHC does not make recommendations regarding radiation protection. It is the customer's responsibility to consult a radiation physicist for advise on radiation protection in x-ray rooms.

The IGS system is equipped with the high performance, highly reliable Performix™ 160A metal X-ray tube which meets requirements for all vascular applications.

Anode diameter	160 mm brazed graphite
Anode rotation	7800 rpm / 130 Hz
Anode Target angle	11.25°
Anode heat storage capacity	2.7 MJ (3.7 MHU)
Anode steady state heat dissipation	6.72 kW
Cathode	Bi-filament design
Coincident focal spot sizes	0.3, 0.6 and 1.0
Fluoroscopic power	<ul style="list-style-type: none"> • 3200 W (continuous) • 4500 W (peak capability for maximum of 10 minutes)
Maximum casing heat storage	5.14 MJ (6.9 MHU)
Continuous casing heat dissipation	3200 W
Maximum anode cooling rate	544 KHU/min (6.72 kW)
Total filtration (IEC 60601-1-3)	1.0 mm Al
Leakage radiation (IEC 60601-1-3)	< 50 mR/h measured at 3.2 kW (125 kV and 25.4 mA)

3.2 Protection against electromagnetic interference hazards

Information below on IEC60601-1-2 and YY9706.102 Electromagnetic Standard Compliance & Documentation can also be found in the IGS System Operator Manual.

General Scope

This equipment complies with IEC60601-1-2 Edition 4.0, IEC60601-1-2 Edition 4.1 and YY9706.102 EMC standard for medical devices.

The IGS system is intended to be used:

- in a PROFESSIONAL HEALTHCARE facility environment,
- in a SPECIAL ENVIRONMENT for IGS systems in OR configuration (vicinity of active HF SURGICAL EQUIPMENT - refer to [Installations Requirements & Environment Control on page 154](#)).

The System is suitable to be used in the electromagnetic environment, as per the limits and recommendations described in the tables here after:

- Emission Compliance level and limits (see [Electromagnetic Emission on page 149](#)).
- Immunity Compliance level and recommendations to maintain equipment clinical utility (refer to [Table 3-2 on page 150](#), [Table 3-4 on page 151](#) and [Table 3-7 on page 154](#)).

Electromagnetic Emission

The IGS System is intended for use in the electromagnetic environment specified below.

The customer or the user of the IGS system should assure that it is used in such an environment.

Table 3-1

Emissions	Compliance	Electromagnetic Environment
Radio-Frequency Emissions CISPR 11 (GB 4824)	Group1 Class A Limits (see the Note below the table).	The IGS system uses Radio Frequency energy only for its internal function. Therefore, its Radio Frequency emissions are very low and are not likely to cause any interference in nearby electronic equipment.
		The IGS system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2 (GB 17625.1)	Not Applicable	The IGS system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions IEC61000-3-3 (GB 17625.2)	Not Applicable	The IGS system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.



NOTE

The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 (GB 4824) class A). If it is used in a residential environment (for which CISPR 11 (GB 4824) class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Electromagnetic Immunity

Electromagnetic Immunity IEC 60601-1-2 and YY9706.102

The IGS system is intended for use in the electromagnetic environment specified below.

The customer or the user of the IGS system should assure that it is used in such an environment.

Table 3-2

IMMUNITY Test	YY9706.102 Test Level	IEC60601-1-2 Ed 4.0 and IEC60601-1-2 Ed 4.1 Test Level (professional health-care environment)	Compliance Level	Electromagnetic Environment
Electrostatic discharge (ESD) IEC61000-4-2 (GB/T 17626.2)	+/- 6 kV contact +/- 8 kV air	+/- 8 kV contact +/- 15 kV air	+/- 8 kV contact +/- 15 kV air	Floors are wood, concrete or ceramic tile or floors are covered with synthetic material and the relative humidity is at least 20 %.
Electrical fast transient/burst IEC61000-4-4 (GB/T 17626.4)	+/- 2 kV for power supply lines +/- 1 kV for input/ output lines 5 kHz burst repetition frequency	+/- 2 kV for power supply lines +/- 1 kV for input/ output lines 100 kHz burst repetition frequency	+/- 2 kV for power supply lines +/- 1 kV for input/ output lines 5kHz & 100 kHz burst repetition frequency	Mains power quality is that of a typical commercial or hospital environment.
Surge IEC61000-4-5 (GB/T 17626.5)	+/- 1 kV line(s) to line(s) +/- 2 kV line(s) to earth	+/- 1 kV line(s) to line(s) +/- 2 kV line(s) to earth	+/- 1 kV line(s) to line(s) +/- 2 kV line(s) to earth	Mains power quality is that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11 (GB/T 17626.11)	<5% U_T (>95% dip in U_T) for 5 sec	0% U_T ; 250/300 cycle	<5% U_T (>95% dip in U_T) for 5 sec 0% U_T ; 250/300 cycle	Mains power quality is that of a typical commercial or hospital environment. If the user of the IGS system requires continued operation during power mains interruptions, it is recommended that the IGS system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC61000-4-8 (GB/T 17626.8)	3 A/m	30 A/m	30 A/m	Power frequency magnetic fields is at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the AC mains voltage prior to application of the test level. 250/300 cycle means 250 periods at 50Hz or 300 periods at 60Hz.				


Table 3-3

IMMUNITY Test	YY9706.102 Test Level	IEC60601-1-2 Ed 4.1 Test Level (professional health-care environment)	Compliance Level	Electromagnetic Environment
Proximity magnetic Fields IEC61000-4-39	-	65 A/m at 134.2 kHz 7.5 A/m at 13.56 MHz	65 A/m at 134.2 kHz 7.5 A/m at 13.56 MHz	It is recommended that magnetic field sources such as RFID readers shall be kept at least 0.15 m from the IGS system.

The IGS system is intended for use in the electromagnetic environment specified below.

The customer or the user of the IGS system should assure that it is used in such an environment.

Table 3-4

IMMUNITY Test	YY9706.102 Test Level	IEC60601-1-2 Ed 4.0 and IEC60601-1-2 Ed 4.1 Test Level (professional health-care environment)	Compliance Level	Electromagnetic Environment
Conducted Radio Frequency IEC61000-4-6 (GB/T 17626.6)	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands ⁽¹⁾	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands ⁽¹⁾	Portable and mobile RF communications equipment is used no closer to any part of the IGS system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated Radio Frequency IEC61000-4-3 (GB/T 17626.3)	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.7 GHz	E1 = 3 V/m ⁽⁴⁾	<p>Recommended separation distance:</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey⁽²⁾, are less than the compliance level in each frequency range.⁽³⁾</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 



NOTE

⁽¹⁾ The ISM (Industrial, Scientific and Medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz; and 40.66 MHz to 40.70 MHz.

**NOTE**

(2) Field strengths from fixed transmitters, such as base stations for cellular telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be performed. If the measured field strength exceeds the RF compliance level above, observe the IGS system to verify normal operation in each use location. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the IGS system.

**NOTE**

(3) Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

**NOTE**

(4) Refer to the table and warning below.

NOTICE

The IGS System is a Large, Permanently-Installed Medical Equipment for which the simulated operation in an anechoic chamber is not feasible and consequently is exempt from the testing requirement specified by IEC 61000-4-3 (GB/T 17626.3).

The IGS System has not been tested for radiated RF immunity over the entire frequency range 80 MHz to 6 GHz.

The IGS System has been tested for radiated RF immunity only at selected frequencies. Use nearby of emitters at other frequencies could result in improper operation.

Table 3-5 YY9706.102 and IEC60601-1-2 field level and frequencies

Tested frequencies (MHz)	Field level (V/m)	Modulation
433.92 (ISM) ⁽⁵⁾	3	80 % AM at 1 kHz rate
915 (ISM) ⁽⁵⁾		
1440		
1750		
1920		
2450 (ISM) ⁽⁵⁾		

**NOTE**

(5) Industrial, Scientific and Medical (ISM) radio bands.

**NOTE**

These are guidelines. Actual conditions may vary.

The associated recommended separation distances as per YY9706.102-2021 are listed in section [Recommended Separation Distances for Portable and Mobile RF Communications Equipment YY9706.102 on page 154](#).

Additional IEC60601-1-2 Ed 4.0 and IEC60601-1-2 Ed 4.1 field level & frequencies - immunity to proximity fields from RF wireless equipment:

Table 3-6 IEC 60601-1-2 Ed4.0 and IEC 60601-1-2 Ed4.1 field level & frequencies

Tested frequencies (MHz)	Field level (V/m)	Modulation
385	27	Pulse modulation (50% duty cycle) – 18 Hz
450	28	Pulse modulation (50% duty cycle) – 18 Hz
710	9	Pulse modulation (50% duty cycle) – 217 Hz
745	9	Pulse modulation (50% duty cycle) – 217 Hz
780	9	Pulse modulation (50% duty cycle) – 217 Hz
810	28	Pulse modulation (50% duty cycle) – 18 Hz
870	28	Pulse modulation (50% duty cycle) – 18 Hz
930	28	Pulse modulation (50% duty cycle) – 18 Hz
1720	28	Pulse modulation (50% duty cycle) – 217 Hz
1845	28	Pulse modulation (50% duty cycle) – 217 Hz
1970	28	Pulse modulation (50% duty cycle) – 217 Hz
2450 (ISM) ⁽⁶⁾	28	Pulse modulation (50% duty cycle) – 217 Hz
5240	9	Pulse modulation (50% duty cycle) – 217 Hz
5500	9	Pulse modulation (50% duty cycle) – 217 Hz
5785	9	Pulse modulation (50% duty cycle) – 217 Hz
5800 (ISM) ⁽⁶⁾	9	Pulse modulation (50% duty cycle) – 217 Hz



NOTE

⁽⁶⁾ Industrial, Scientific and Medical (ISM) radio bands.



NOTE

These are guidelines. Actual conditions may vary.

Equipment used for tests:

- RF signal generator.
- RF power amplifier.
- Transmitting antenna.
- Field sensor.
- Field meter.

WARNING



PORTABLE RF COMMUNICATIONS EQUIPMENT INCLUDING PERIPHERALS (SUCH AS ANTENNA CABLES AND EXTERNAL ANTENNAS) SHOULD BE USED NO CLOSER THAN 30 CM (12 INCHES) TO ANY PART OF THE IGS SYSTEM INCLUDING CABLES SPECIFIED BY THE MANUFACTURER. OTHERWISE, DEGRADATION OF THE PERFORMANCE OF THIS EQUIPMENT COULD RESULT.

Recommended Separation Distances for Portable and Mobile RF Communications Equipment YY9706.102

Table 3-7

Frequency of Transmitter	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
Equation	$d=[3.5 / \sqrt{P}] \cdot \sqrt{P}$	$d=[3.5 / \sqrt{E1}] \cdot \sqrt{P}$	$d=[7 / \sqrt{E1}] \cdot \sqrt{P}$
Rated Power of Transmitter (watts)	Distance (meters)	Distance (meters)	Distance (meters)
10 mW	0.11	0.11	0.22
100 mW	0.37	0.37	0.74
1	1.1	1.1	2.3
10	3.7	3.7	7.4
100	12	12	23

For transmitters rated at a power not listed above, the DISTANCE can be estimated using the equation in the corresponding column, where P is the power rating of the transmitter in watts (W) according to the transmitter manufacturer.



NOTE

These are guidelines. Actual conditions may vary.

Limitations Management

Adhering to the distance separation recommended in section [Recommended Separation Distances for Portable and Mobile RF Communications Equipment YY9706.102 on page 154](#), between 150 kHz and 2.5 GHz, will reduce disturbances recorded at the image level, but may not eliminate all disturbances. However, when installed and operated as specified herein, the system will maintain its essential performance by continuing to acquire, display, and store diagnostic quality images safely.

For example, a 1 W mobile phone (800 MHz to 2.5 GHz carrier frequency) shall be put 2.3 meters (refer to [Table 3-7 on page 154](#) in section [Recommended Separation Distances for Portable and Mobile RF Communications Equipment YY9706.102 on page 154](#), apart from the IGS system (in order to avoid images interferences risks).

Installations Requirements & Environment Control

WARNING



USE OF ACCESSORIES, TRANSDUCERS AND CABLES OTHER THAN THOSE SPECIFIED OR PROVIDED BY THE MANUFACTURER OF THIS EQUIPMENT COULD RESULT IN INCREASED ELECTROMAGNETIC EMISSIONS OR DECREASED ELECTROMAGNETIC IMMUNITY OF THIS EQUIPMENT AND RESULT IN IMPROPER OPERATION.

In order to minimize interference risks, the following requirements shall apply:

- Separated Power supply distribution panel & separated power line.
- Stacked components and equipment:
The IGS system should not be used adjacent to or stacked with other equipment.
- Low frequency magnetic field: other electrical equipment may disturb and interfere with these IGS components. The control of the clearing distances from the noise sources is recommended from the HF electrosurgery generator, power supplies converters from nearby monitors or from other close electrical equipment). Refer to respective device manufacturers instructions and recommendations in such cases.
- In order to minimize the risk of interference from a nearby MRI device, the recommended maximum static magnetic field amplitudes are specified below:
 - Static field is specified less than <1 Gauss in Exam room, and in the Control Room.
 - Static field is specified less than <3 Gauss in the Technical Room.
- Electrostatic discharges environment and recommendations:
 - In order to reduce electrostatic discharge interference, install a charge dissipative floor material to avoid electrostatic charge buildup.
 - The relative humidity shall be at least 20 percent as defined in [4.1.1 Humidity on page 156](#).
 - The dissipative material shall be connected to the system ground reference, if applicable.

4 Environmental Requirements

4.1 Humidity, Temperature and Altitude

4.1.1 Humidity

Table 4-1 Relative Humidity (non- condensing)

	MIN	MAX
Exam Room	20%	70%
Control Room	20%	75%
Technical Room	20%	75%

4.1.2 Temperature and Altitude

The system is certified for use up to 3000 m. The permissible atmospheric pressure conditions of use are between 700 hPa and 1060 hPa.

Above 2000 m, the thermal dissipation is reduced because the air pressure is lower. Therefore, a temperature derating shall be applied for the Technical Room as defined in the table below.

The CE Fluoro UPS can be used only below 2000 m.

Table 4-2 Exam Room and Control Room - Temperature

	MIN	MAX	RECOMMENDED
Exam Room	+15°C (+59°F)	+32°C (+90°F)	+20°C (+68°F)
Control Room	+15°C (+59°F)	+35°C (+95°F)	+20°C (+68°F)

Table 4-3 Technical Room - Temperature

	Temperature up to 2000 m			Temperature above 2000 m		
	MIN	MAX	RECOMMENDED	MIN	MAX	RECOMMENDED
Technical Room (with 1 kVA)	+10°C (+50°F)	+30°C (+86°F)	+20°C (+68°F)	+10°C (+50°F)	+20°C (+68°F)	+20°C (+68°F)
Technical Room (with 8 kVA or the Fluoro UPS)	+15°C (+59°F)	+25°C (+77°F)	+20°C (+68°F)	+15°C (+59°F)	+20°C (+68°F)	+20°C (+68°F)

**NOTE**

For the systems that are planned to be installed at the second floor or above, the temperature and humidity of the rooms that are directly below the gantry room should be the same as the Exam Room requirement.

Differences in temperature or humidity between the Exam room and the room located below will cause condensation within the gantry or patient table, resulting in part failure or rust. Failure to do so will void the equipment warranty. Avoid above grade installations if the temperature is high in the area below the cables entrance of the gantry or table.

4.2 Heat Output

In the table:

- Moderate Use corresponds to 8 cases per a 10 hours day,
- Typical Use corresponds to 11 cases per a 10 hours day,
- Maximum Use is maximum peak power during exam.

Table 4-4

		HEAT OUTPUT							
		Stand by		Moderate Use		Typical Use		Maximum Use	
Room	Core System	kW	BTU/hr	kW	BTU/hr	kW	BTU/hr	kW	BTU/hr
Exam Room	Gantry and Table	0.41	1,399	0.55	1,877	0.89	3,037	1.62	5,528
	6 19" monitors on suspension or	0.30	1,024	0.30	1,024	0.30	1,024	0.30	1,024
	LDM suspension with 2 backups	0.10	341	0.10	341	0.10	341	0.10	341
	Typical Injector	0.09	307	0.09	307	0.09	307	0.09	307
Control Room	DL console and Live monitor	0.10	341	0.10	341	0.10	341	0.10	341
Technical Room	C-FRT Cabinet	0.70	2,388	1.02	3,480	1.53	5,221	2.16	7,370
	PDU	0.50	1,706	0.50	1,706	0.50	1,706	0.50	1,706
	Tube Chiller	2.53	8,633	4.49	15,321	5.49	18,733	6.93	23,646
	Detector Conditioner	0.21	717	0.21	717	0.21	717	0.21	717
	UPS 1 kVA	0.15	500	0.15	500	0.15	500	0.15	500
	UPS 8 kVA	0.52	1,774	0.52	1,774	0.52	1,774	0.52	1,774
	Fluoro UPS UL	2.14	7,302	2.14	7,302	2.14	7,302	2.14	7,302
	Fluoro UPS CE	0.585	2,000	0.585	2,000	0.585	2,000	0.585	2,000
Total for Core System with the 1 kVA UPS		4.79	16,332	7.21	24,590	9.06	30,902	11.86	40,456
Total for Core System with the 8 kVA UPS		5.16	17,606	7.58	25,864	9.43	32,176	12.23	41,730
Total for Core System with the Fluoro UPS UL		6.78	23,134	9.20	31,392	11.05	37,704	13.85	47,258
Total for Core System with the Fluoro UPS CE		5.225	17,832	7.645	26,090	9.495	32,402	12.295	41,956

4.3 Acoustic Specifications

- Less than 50 dB (A) at 1 meter for Gantry.
- Limited to 50 dB (A) at 1 meter for Omega IV Table and Omega V Table.
- Limited to 58 dB (A) at 1 meter for Innova^{IQ} OR Table.
- Limited to 59 dB (A) at 1 meter for C-FRT Cabinet and NPA PDU.
- Limited to 60 dB (A) at 1 meter for the Tube Chiller.
- Limited to 52 dB (A) (background of 35 dB (A)) at 1 meter for Detector Conditioner.
- Limited to 39 dB (A) at 1 meter for UPS 8 kVA Gen1.
- Less than 37 dB (A) at 1 meter for UPS 8 kVA Gen2.
- Less than 60 dB (A) at 1 meter for the Fluoro UPS.
- less than 55 dB (A) at 20 degrees Celsius, measured in the operators head position, 20 cm in front of the keyboard's right corner, at 1.30 m above the floor, and in a distance of 1 meter at all four sides.

4.4 Room Light

4.4.1 Requirements for Lighting

Requirement for lighting concern the following, general, light-technique characteristics:

- Illuminator level.
- Lighting distribution.
- Preventing the operator from being dazzled by the light (by direct light sources or by reflection on bright objects).

The Illumination level must be compliant with established lighting technical rules and be as constant as possible.

Technical Room, Exam Room and Control Room shall be provided with appropriate lighting in the maintenance area (maintenance area to be considered are service workplaces). It corresponds to service areas as defined for any of the product components.

The minimum required average luminance E_m shall be of 500 lux and minimum color rendering factor R_a of 80 as per IEC/EN 12464-1 (Light and lighting. Lighting of work places. Indoor work places: Illumination requirements for indoor workplaces corresponding to assembly of medium size electrical components, e.g. control panel) for the electrical industry).

4.4.2 Windows and Curtains

When the Exam Room has a window with an aperture outside of the controlled light area (day light, other...) a curtain has to maintain the light intensity under a limit fixed to 150 lux.



NOTE

In Germany: Ambient luminance of 100 lux maximum is required to maintain Exam Room class 2 according to DIN 6868-157.

4.4.3 Surgical Lights

WARNING

If a surgical light is installed by the customer, it has to be powered from an independent power supply (provided by the hospital not by the System).

5 Electrical Requirements

5.1 System Electrical Ratings

Table 5-1

Nominal voltage	Frequency	Type of power input		Power consumption				
		With 8 kVA UPS	With 20 kVA UPS	System Off	Standby	Long time	Momentary	Peak
380 V	50 Hz or 60 Hz	3~	3N~	0.4 kVA	3.1 kVA	15 kVA	100 kVA	150 kVA
400 V								
415 V								
440 V	60 Hz		Not compatible					
480V			3N~					

Long time rating is measured in fluoroscopy mode at 30 fps, 120 kV, 89 mA, 10 ms.

Momentary rating is measured in record DSA mode at 7.5 fps, 125 kV, 640 mA, 50 ms.


For the rating of the external devices not powered by the system (AW, injector, and so on), refer to the OEM documentation.

5.2 Power Distribution Schematics

Information below specifies the cables provided by GE and the cables provided by the Hospital. Refer to MDP [5.3 Cabling Requirements on page 164](#).

5.2.1 System with 1 or 8 kVA UPS

Figure 5-1 Power Distribution with 1 or 8 kVA UPS




NOTE

UPS 1 kVA is only provided with System with Omega Table.

5.2.2 System with Fluoro UPS

Figure 5-2 Power Distribution with Fluoro UPS



5.2.3 System with Fluoro UPS and IT Electrical Network

The Fluoro UPS requires a Neutral line connected to the Protective Earth. For hospitals with an IT Electrical Network, a transformer is required with Delta-Wye or Delta-Star connection.


This transformer shall be provided by the customer; its characteristics shall be:

- 30 kVA minimum.
- Secondary star 3 Ph+N.
- The power distribution shall be of TNS type with the Neutral grounded.
- The transformer impedance shall be 4.5% or less (this parameter is also called %Z or short circuit voltage).

It is also the responsibility of the customer to provide:

- The box of the transformer to avoid access to live parts according to local regulations.
- The current protections at the output of the transformer, fuses or breaker, as per local regulations. The suggested rating is 50A.

Figure 5-3 Power Distribution with Fluoro UPS and IT Electrical Network



5.2.4 Power distribution for 3rd party monitors on GE 19" monitors suspensions

GE 19" monitors suspensions are pre-equipped for the connection of one monitor for 3rd party devices (ECG, ultrasound,...). It allows the connection of the video signal through a VGA cable or a DVI extender. The power connection is pre-routed between the connection box in the suspension and the monitor.

This 3rd party monitor shall not be powered by the system, but by a dedicated hospital power outlet. This power outlet shall be fitted with a dedicated overcurrent protection (e.g. circuit breaker) and a means of isolation with provision for Lockout/Tagout, accessible at all times. A fuse or a semiconductor device shall not be used as a means of isolation. A unique device can serve as circuit breaker and means of isolation, provided it meets all the requirements below.

The installation shall be done in accordance with all local regulations.

1. Technical requirements for the overcurrent protection:

- The rating of the circuit breaker shall be less than 12A.
- The breaking capacity of circuit breaker must be adapted to upstream input line breaking capacity.

2. Technical requirements for the means of isolation:

The means of isolation shall comply with the following requirements:

- rated Impulse Withstand Voltage (Uimp): 4 kV
- provision for Lockout/Tagout (LOTO)
- opens all poles simultaneously
- direction of movement of the actuator complies with IEC 60447 (The OFF to ON direction must be left to right, bottom to top or clockwise)

The power cable between the hospital outlet and the connection box of the suspension shall be provided by the hospital, its diameter shall be in accordance with the rating of the overcurrent protection.

5.3 Cabling Requirements

It is the customer's responsibility to ensure that the electrical installation is compliant with local regulations, such as NFPA99 (Health Care Facilities Code) or 60364-7-710 (Requirements for special installations or locations - Medical locations).

To avoid risk of electric shock, this equipment must only be connected to a mains power supply with Protective Earth.

The power supply and ground cables shall be dedicated to the system. They must not be used to supply other systems. Power supply and ground cables shall be kept separated from other room System cables and must be connected to the same distribution panel. They must run near one to the other.

The power cables, ground cables and EPO cables provided by the customer shall be compliant with local regulations (e.g. UL, NFPA 70, CSA, IEC, CCC).

5.3.1 Power Source Characteristics

5.3.1.1 Transformer

The transformer powering the system shall meet the following requirements:


- 150 kVA minimum for input voltage of 380 V and 400 V.
- 100 kVA minimum for input voltage of 415 V and 480 V.
- The transformer impedance shall be 4.5 % or less (this parameter is also called %Z or short circuit voltage).

The protection of the output of the transformer shall be sized 100A minimum

It shall provide a neutral for the systems with the 20kVA UPS (else refer to [5.2.3 System with Fluoro UPS and IT Electrical Network on page 162](#)).

5.3.1.2 Power Cables (Feeders)

Figure 5-4



The power cables shall meet all the requirements below, whichever is the most stringent:

- 2 AWG / 35mm² minimum.
- The size of the power cables shall be adapted to the protection of the transformer, this means that the ampacity of the cables shall be greater or equal to the rating of the protection.
- The DC resistance of the power cables shall be less than the values defined in the [Table 5-2 on page 165](#). This means that the DC resistance of the cable 1a + cable 1b + cable 1c, cable 2a + cable 2b + cable 2c and cable 3a + cable 3b + cable 3c in [Figure 5-4 on page 165](#) shall all be less than the “Per phase” value.

The protective Earth cables shall not be smaller than the power cables (cables #1 & #2). The Protective Earth cable between the MDP and the PDU (cable #2) shall be separated from the power cable.

Max Line Impedance at the entry of the X-rays Generator in C-FRT Cabinet:

Table 5-2

Mains Voltage (V)	380	400	415	480
Max Line Impedance Phase to phase (Ω)	0.09	0.096	0.102	0.12
Max Line Impedance Per phase (Ω)	0.045	0.048	0.051	0.06

The maximum cable length between the MDP and the PDU shall be 20m (65.6ft). This cable shall be a copper cable and the cable insulation temperature shall be 90°C. Use cable Type S for North America



NOTE

The cables 1c, 2c and 3c are provided with the system, and are AWG2, 8.2m long. The DC resistance of each cable is 0.0042 Ohms.

The [Table 5-3 on page 165](#) is provided for information for the calculation of the DC resistance of copper cables, the calculation shall be done using the actual values in the datasheet of the cables that will be used.

Table 5-3

Resistivity of copper wires (for information)		
AWG	Diameter (mm ²)	Resistivity (Ω/km)
0000 (4/0)	107	0.1608
000 (3/0)	85	0.2028

Table 5-3 (Table continued)

Resistivity of copper wires (for information)		
AWG	Diameter (mm ²)	Resistivity (Ω/km)
00 (2/0)	67	0.2557
0 (1/0)	53	0.3224
1	42	0.4066
2	35	0.5127

5.3.2 Fluoro UPS Cables

5.3.2.1 Fluoro UPS UL Power Cables

The minimum gauge of the power and Protective Earth from the MDP to the FUPS and from the FUPS to the PDU shall be 10 mm² / 6 AWG (cables # 3 and # 4 on [Figure 5-2 Power Distribution with Fluoro UPS on page 162](#)).

The Protective Earth cables shall not be smaller than the power cables.

5.3.2.2 Fluoro UPS CE Power Cables

The gauge of the power and Protective Earth from the MDP to the FUPS shall be 10 – 16 mm² / 7 AWG - 5 AWG (cables # 3 on [Figure 5-2 Power Distribution with Fluoro UPS on page 162](#)).


The gauge of the power and Protective Earth from the FUPS to the PDU shall be 10 mm² / 7 AWG. (cable # 4 on [Figure 5-2 Power Distribution with Fluoro UPS on page 162](#)).

The Protective Earth cables shall not be smaller than the power cables.



NOTE

The entry point of the power cables at the back of the UPS is located 85 cm (2.8 feet) from the floor.

Figure 5-5

Dimension in cm (feet)

5.3.2.3 Fluoro UPS Ethernet Cable

A 6 m Ethernet cable between the C-FRT Cabinet and the Fluoro UPS is provided with the system. If a longer cable is needed, it shall be provided by the hospital; it shall be Cat5 minimum.

5.3.3 EPO Cables and additional buttons

5.3.3.1 EPO Cable between the MDP and the PDU

The hospital shall provide an EPO cable between the MDP and the PDU; its minimum gauge shall be 1 mm² / 16 AWG and shall be in accordance with the rating of the fuse F2 of the MDP.

A 12 m EPO cable between the PDU and the Fluoro UPS is provided with the system. If a longer cable is needed, it shall be provided by the hospital, its minimum gauge shall be 1 mm² / 16 AWG.

5.3.3.2 Additional EPO buttons

The PDU is provided with an EPO button on its front panel and provides the connection for additional EPO buttons (in Exam Room or Control Room).

The customer is responsible for the procurement, delivery and installation of the cables and EPO buttons. It is recommended not to locate the EPO in a high traffic area, or an area where someone could accidentally bump it, press it, and mistake it for a different button.

The EPO buttons shall provide 2 Normally Closed contacts, compatible with 24 V AC and in accordance with the MDP transformer rating. The maximum length of the cables shall be 24 m, the recommended diameter is 1 mm² / 16 AWG.

Once activated, the EPO button shall require a user action to deactivate it (for instance "Push to activate - Push to release" or lever type).

WARNING



ACCIDENTAL ACTIVATION OF AN EPO BUTTON CAN LEAD TO UNINTENDED INTERRUPTION OF A PATIENT PROCEDURE COUPLED WITH A DELAY WHEN RESETTING THE EPO, THEN REPOWERING AND REBOOTING THE SYSTEM.

THE POSITION OF THE ADDITIONAL EPO BUTTONS IN THE EXAM ROOM AND CONTROL ROOM SHALL ALLOW TO MINIMIZE ACCIDENTAL ACTIVATION BY A USER, SHOCK WITH A GURNEY...

PENDING HOSPITAL POLICIES AND LOCAL REGULATIONS, IT IS RECOMMENDED THAT THE EPO BUTTON BE PROTECTED AGAINST ACCIDENTAL ACTIVATION (FOR INSTANCE: RECESSING THE EPO INTO A WALL OR PUTTING IT INTO A BOX/FRAME (FIG BELOW) AND LABELED TO DISTINGUISH IT FROM OTHER BUTTONS, FOR INSTANCE: ROOM LIGHT BUTTON, CODE BLUE BUTTON, etc.) WITHIN THE ROOM.



Legend of Figure 5-6 EPO Schematic with 1 kVA or 8 kVA UPS on page 169 and Figure 5-7 EPO Schematic with Fluoro UPS 20 kVA on page 169:

Item	Description
	Black line: Internal wiring and GE responsibility wiring
	Red line: External wiring (Customer responsibility)
[1]	Remote EPOs
[2]	UPS EPO
[3]	PDU
[4]	PDU EPO
[5]	PDU Transformer temperature sensor
[6]	MDP

Item	Description
[7]	MDP EPO
[8]	MDP EPO internal circuit

Figure 5-6 EPO Schematic with 1 kVA or 8 kVA UPS





Figure 5-7 EPO Schematic with Fluoro UPS 20 kVA



**NOTE**

J35 connection:

- UPS 1 kVA: pair 1, 2 connected to Transformer EPO output, pair 3, 4 is shorted.
- UPS 8 kVA: pair 1, 2 connected to Transformer EPO output, pair 3, 4 is shorted.
- UPS 20 kVA (Fluoro UPS): pair 3, 4 connected to Transformer EPO output, pair 1, 2 is shorted.

5.4 Mains Disconnect Panel

5.4.1 General Information

5.4.1.1 Introduction

The Mains Disconnect Panel (MDP) is the electric panel which is the interface between the Hospital mains and the System. It allows the power connections from the hospital power to the input of the PDU of the system and to the Fluoro UPS if present. It provides the LOTO (lock out – tag out) functions that allows safe service operation, and is part of the EPO (Emergency Power Off) function.

As the requirements applicable to electric panels vary from a country to another, information below lists the GE mandatory requirements to provide safe system operation and the installation precautions, in addition to the local regulatory requirements.

Information given shall allow the Customer to build the MDP in compliance with GE's rules. In addition, the following MDP can be ordered through the GE accessory catalog:

- MDP CE (E46001BB), certified IEC 61439-1,
- MDP UL (E46001BA), certified UL508A for USA,
- MDP UL (E46001BD), certified UL508A for USA and OSPHD.



The Customer MDP is not covered by the GEHC product certification. The association of the Allia™ System and the Customer MDP is not covered by the GEHC product certification.

GE specifically disclaims any and all liability arising out of or relating to the use or performance of the MDP and the cables in the scope of the System Pre-Installation Manual, including, and without limitation, any liability or claims relating to patient injury, death, or the reliability of such MDP.

The mechanical and electrical installation of the MDP is fully under the customer and the installer responsibility.

The customer is responsible for ensuring that all requirements from the System Pre-Installation Manual are met.

5.4.1.2 Pre-Installation

It is the customer responsibility to ensure that the MDP and its input and output cables are installed prior to the GE equipment (PDU, other cabinets, FUPS option, etc.) to ensure that standard GE Service Process can be followed during the System installation. The connection of the MDP to the GE equipment shall only be made in presence of a GE Service representative.

It is recommended that the vendor contacts GE Service representative and reviews the site planning details before the MDP is installed.

**NOTE**

GE will not be responsible for any delay in installation if the MDP is not mounted and its cables not routed before GE parts arrive on site.

5.4.1.3 Spare Parts

The customer is responsible for providing and replacing any part of MDP.

5.4.2 Mandatory Construction Requirements

5.4.2.1 Input Power

The MDP shall be functional within one the following input voltage and frequency ranges from the Hospital mains:

- Voltage range for systems without the Fluoro UPS:
 - 380 V +/-10% 3~, 50 Hz or 60 Hz +/-3 Hz
 - 400 V +/-10% 3~, 50 Hz or 60 Hz +/-3 Hz
 - 415 V +/-10% 3~, 50 Hz or 60 Hz +/-3 Hz
 - 440 V +/-10% 3~, 60 Hz +/-3 Hz
 - 480 V +/-10% 3~, 60 Hz +/-3 Hz
- Voltage range for systems with the Fluoro UPS:
 - 380 V +/-10% 3N~, 50 Hz or 60 Hz +/-3 Hz
 - 400 V +/-10% 3N~, 50 Hz or 60 Hz +/-3 Hz
 - 415 V +/-10% 3N~, 50 Hz or 60 Hz +/-3 Hz
 - 480 V +/-10% 3N~, 60 Hz +/-3 Hz

5.4.2.2 Breakers

The MDP shall provide a main breaker at its input, its specifications shall be:

- Current Rating: 100 A.
- It shall be capable of withstanding an inrush current of 2000 A for 10 ms.
- The voltage rating shall be the MDP nominal input line voltage +10%: i.e., 380 V + 10%, 400 V + 10%, 415 V + 10%, 440 V + 10% or 480 V +10%.
- The frequency range shall be adapted to the input line frequency i.e., 50 Hz +/-3 Hz or 60 Hz +/-3 Hz.
- The Short Circuit Current Rating (SCCR) shall be adapted to the input line source short circuit capacity.

This command of this breaker shall be accessible from the outside of the MDP, in order to be able to rearm it without opening the MDP after an emergency power off.

For systems with the FUPS, the MDP shall provide a second breaker for protection of the FUPS input. Its specifications shall be:

- 3 poles type
- Voltage rating: same as the MDP main breaker
- Frequency range: same as the MDP main breaker

- Current Rating: 50 A
- Short Circuit Current Rating (SCCR): 20 kA

For systems with the Fluoro UPS, the breaker for the FUPS input protection shall be powered by the MDP main breaker.

5.4.2.3 Terminal Blocks

The MDP shall have an input mains terminal block rated in accordance with the hospital input voltage. It shall be capable of holding minimum 35 mm² / 2 AWG cable for the 3 phases, protective Earth and neutral (only for systems with the Fluoro UPS).

The MDP shall provide an output terminal block rated in accordance with the MDP input voltage to connect the output from the MDP main breaker to the system. This terminal block shall be capable of holding minimum 35 mm² / 2 AWG cable for the 3 phases.

For systems with the Fluoro UPS, the MDP shall have a terminal block to connect an input neutral from the mains and an output neutral to the Fluoro UPS, and it shall have an output terminal block rated to the hospital input voltage to connect the mains input power from the MDP to the FUPS. This terminal block shall be capable of holding minimum 10 mm² / 6 AWG cable for the 3 phases and neutral.

5.4.2.4 Protective Earth

The MDP shall have a ground bar / ground terminal to connect the protective Earth cables:

- from the hospital mains,
- to the system,
- to the FUPS (if present).

5.4.2.5 Indicators

The MDP shall have lights to indicate the presence of voltage. The presence of voltage on each input line shall be indicated by at least having lamps between Line1-Line2 and Line2-Line3. The recommended color for these lamps is green.

5.4.3 Mandatory EPO Requirements

The MDP shall provide an emergency power off (EPO) button on its front.

The EPO button shall not be of momentary type.

The EPO button shall have 2 NC contacts:

- one NC contact is to trip the MDP input breaker,
- the other NC contact is to activate the UPS EPO input (1 kVA or 8 kVA or FUPS) to turn off the UPS output (this connection is done inside the PDU).

The MDP shall provide a terminal block to connect external cables to the 2 NC contacts of the MDP EPO.

When the MDP EPO or the PDU EPO is pressed, the MDP shall not provide any output voltage without any additional action on the EPO buttons and on the MDP input breaker.

The EPO button shall be protected against accidental activation, in order to prevent from accidental power OFF as shown below or equivalent.

Figure 5-8 EPO Button




5.4.4 Mandatory LOTO Requirements

The MDP shall provide means of disconnecting the mains power from the system, with LOTO capability to ensure safe service operation. It can be done by the input breaker if it has disconnecting capability, or by a separate disconnection device.

An operator should be able to apply LOTO without opening the MDP box. When a LOTO device is installed on the MDP input breaker or on the disconnecting device, there shall be no voltage at the output of the MDP.

5.4.5 Mandatory Labeling Requirements

The input mains terminal block and the output terminal blocks of the MDP shall be labeled to indicate the 3 lines as shown below or equivalent:



For systems with the Fluoro UPS, the MDP input neutral terminal block and the neutral output terminal block to the Fluoro UPS shall be labelled with the IEC 60445 symbol as shown below or equivalent:



The ground bar shall be marked with the IEC 60417-5019 symbol as shown below:



5.4.6 Other Mandatory Requirements

The MDP and the external cables shall be compliant to all applicable local regulations, in particular to the standards applicable to Industrial Control Panels or Low-voltage switch gear and control gear assemblies, such as UL508A for USA or IEC 61439-1 for Europe.

The MDP enclosure shall be grounded if its enclosure is metallic, and there shall be no access to hazardous voltages. The enclosure shall provide enough rigidity to avoid hazardous situations in case of shock or impact and shall be designed in accordance with the local regulations.

Local regulation may require the MDP to have a door interlock mechanism to prevent from opening the door when the main breaker is on.

5.4.7.2 Minimum Components Specifications for the CE MDP

Table 5-4


Component	Label (refer to Figure 5-9 CE MDP - Power and Control on page 174)	Rating
Input Circuit Breaker	CB1	4 Pole, 100A, 20 kA SCCR and adapted to the input line, Vin + 10% 50 Hz or 60 Hz 2000 A inrush current withstand capability for at least 10 ms
Circuit Breaker	CB2	50 A, 20 kA SCCR, Vin+ 10%, 50 Hz or 60 Hz
Fuse	F1	2A Time delay, Vin+10% Based on green indicator lights power ratings
Fuse	F2	2A, 24 VAC+10% Based on transformer power rating and transformer load current rating
Fuse	F3	1A Time delay, Vin+10% Based on transformer power rating and transformer input current
Time delay relay	T2	24 VAC+10% Shall have 1 NO contact Time delay setting shall be min 2 Sec
Auxiliary relay	R1, R2*	24 VAC+10% Shall have 1 NC contact, 1 NO contact * R1 and R2 part numbers shall be identical (same manufacturer, same reference) Preferred components for R1 and R2 relay are: <ul style="list-style-type: none"> • GE: PRC1S13-BDL • ABB: CR-M024AC2L • Schneider Electric: 782XBXM4L-24A • Omron: MK2PI
Shunt Trip	ST	24 VAC+10% Shunt trip opens the MDP input main breaker when the shunt trip is energized
2 Pilot lights Green	-	Vin+10%, 50 Hz or 60 Hz
Transformer	TR1	Power rating: 50 VA or based on power ratings of components used at transformer output Input: 380 VAC or 400 VAC or 415 VAC Output: 24V Frequency: 50 or 60 Hz Double insulation as per standard IEC61558 The sum of power ratings of R1, shunt trip and timer shall be less than transformer power rating
EPO	-	Mushroom button with 2 NC contacts Rated for 24 VAC, 50 mA

Table 5-4 (Table continued)

Component	Label (refer to Figure 5-9 CE MDP - Power and Control on page 174)	Rating
Cable for MDP internal Control circuitry	-	Min 1 mm ² and in accordance with the fuses rating


5.4.7.3 Recommended UL Schematics

Figure 5-11 UL MDP - Power and Control



NOTE
 Neutral is not required by Imaging system.
 Neutral is required when Fluoro UPS is used.

Figure 5-12 UL MDP - I/O Interfaces



5.4.7.4 Minimum Components Specifications for the UL MDP

Table 5-5

Component	Label (refer to Figure 5-11 UL MDP - Power and Control on page 176)	Rating
Input Circuit Breaker	CB1	3 Pole, 100 A, 20 kA SCCR and adapted to the input line, 480 VAC+10% 60 Hz 2000 A inrush current withstand capability for at least 10 ms
Circuit Breaker	CB2	50 A, 20 kA SCCR, 480 VAC+10%
Fuse	F1	2A Time delay, 480 VAC+10% Based on green indicator lights power ratings
Fuse	F2	2A, 24 VAC+10% Based on transformer power rating and transformer load current rating
Fuse	F3	1A Time delay, 480 VAC+10% Based on transformer power rating and transformer input current
Time delay relay	T2	24 VAC+10% Shall have 1 NO contact Time delay setting shall be min 2 s
Auxiliary relay	R1, R2*	24 VAC+10% Shall have 1 NC contact, 1 NO contact * R1 and R2 part numbers shall be identical (same manufacturer, same reference) Preferred components for R1 and R2 relay are: <ul style="list-style-type: none"> • GE: PRC1S13-BDL • ABB: CR-M024AC2L • Schneider Electric: 782XBXM4L-24A • Omron: MK2PI
Shunt Trip	ST	24 VAC+10% Shunt trip shall open the MDP input breaker when shunt trip is energized
2 Pilot lights Green	PH1-PH2 PH2-PH3	480 VAC +10%
Transformer	CP T1	Power rating: 50 VA or based on power ratings of components used at transformer output Input: 480 VAC Output: 24 VAC Frequency: 60 Hz +/-3 Hz Double insulation as per UL 5085-1 standard The sum of power ratings of R1, shunt trip and timer shall be less than transformer power rating

Table 5-5 (Table continued)

Component	Label (refer to Figure 5-11 UL MDP - Power and Control on page 176)	Rating
EPO	-	Mushroom button with 2 NC contacts Rated for 24 VAC, 50 mA
Cable for MDP internal Control circuitry	-	Min 16 AWG and in accordance with the fuses rating

5.4.8 Checklist

The following checklist shall be filled and given to the Field Engineer before connecting the MDP to the system.

Table 5-6

Test	Expected Result	OK / NOK
Functional Tests		
Initial state: the MDP main breaker is off, power is available at its input. A jumper is installed between TB2 7 & 8 Turn on the MDP main breaker.	The indicator lights on MDP front panel are ON.	
	The voltage at TB3 is the same as the MDP input voltage.	
	For systems with the FUPS, the voltage at TB4 is the same as the MDP input voltage.	
Press the EPO push button on MDP front panel.	The indicator lights on the MDP front panel are turned off.	
	The MDP main breaker is opened.	
	There is no voltage at TB3.	
	For systems with the FUPS, there is no voltage at TB4.	
The dry contact between TB2 5 & 6 is open.		
Check it is possible to apply the LOTO on the MDP input breaker or on the disconnecting device.	It is possible to apply the LOTO on the MDP input breaker or on the disconnecting device.	
Documentation		
Check a LOTO procedure is provided with the MDP.	The LOTO procedure is present.	
Components Ratings		
Check that the components ratings are compliant with the requirements of 5.4.7.2 Minimum Components Specifications for the CE MDP on page 175 or 5.4.7.4 Minimum Components Specifications for the UL MDP on page 177 .	The components ratings are compliant with the requirements of 5.4.7.2 Minimum Components Specifications for the CE MDP on page 175 or 5.4.7.4 Minimum Components Specifications for the UL MDP on page 177 .	

5.5 External Interfaces


5.5.1 X-Ray ON lights

NOTICE

The X-Ray ON lamp must be installed in the Exam Room in conformity to the standard IEC/EN 60601-2-43. The X-Ray ON lamp shall be visible by the operator in all the locations defined for the personnel who may receive scattered radiation.

The system provides a dry contact to trigger a low voltage relay (24 V max) that drives the X-Ray ON lights. The customer is responsible for the procurement, delivery and installation of the power supplies, relay, cables and the X-Ray ON lights.

Figure 5-13




The cables are connected to the PDU on an open contact. The diameter of the cables shall be 2 mm² maximum.

5.5.2 System ON light

The system provides a dry contact to trigger a low voltage relay (24 V max) that can drive a System ON light. The customer is responsible for the procurement, delivery and installation of the power supplies, relay, cables and the System ON light.

Figure 5-14




The cables are connected to the PDU on an open contact. The diameter of the cables shall be 2 mm² maximum.

5.5.3 Room lights

The system provides a dry contact to trigger a low voltage relay (24 V max) that can drive the Exam Room lights. The customer is responsible for the procurement, delivery and installation of the power supplies, relay, cables and the room lights.

Figure 5-15




The cables are connected to the PDU on a closed contact. The diameter of the cables shall be 2 mm² maximum.

5.5.4 Room door interlock

The system provides a room door interlock that can prevent X-Ray emission when the door is open. The IEC 60601-2-43 requires not to install door interlocks. It is the responsibility of the installer to verify that the connection of this interlock is not in contradiction with local regulation. In case of conflict, the local regulation shall prevail.

This switch shall be closed when the door is closed, it shall be compatible with 24 V DC.

Figure 5-16




To disable the door interlock: the pins 2 and 3 from J436 shall be shorted. The diameter of the cables connected to the cabinet shall be 2 mm² maximum.

5.5.5 Video distribution

The system can provide DVI outputs (1280 x 1024, 60 Hz) for each of its 3 displays (live, roadmap and review). Only 4 streams of each display can be provided (including the images displayed on the LDM). These video links are 36 meters optical cables. The power supply (PSU) of the optical receiver (Rx) shall not be powered by the system.

In case more than 4 video streams of a display are needed, additional DVI video splitters and DVI cables shall be provided and installed by the customer as on schematic below. These splitters shall not be installed inside the cabinets, nor be powered by the system.

Figure 5-17



With the LDM option, a 2MP copy of the LDM image (DVI 1920 x 1080, 60 Hz) can be provided as an option, through a 36 m optical cable. The monitor shall be provided and installed by the customer. The monitor and the optical receiver shall not be powered by the system.

5.6 Additional Full UPS

If it is required to power continuously the system in record mode during power failure, a 150 kVA UPS can be used in front of the system. Such an UPS will provide to the customer about 10 minutes of autonomy. This UPS comes in addition to the UPS provided with the system.

In case the system is powered by a mains voltage of 380, 400 or 415Vac (50 or 60 Hz) it is also possible to use the UPS 93E-G2 120 kVA from Eaton (E45011ND) for this purpose.


5.7 Grounding Schematics



NOTE


For Third party monitor suspension grounding, refer to 5.8 Third-Party Monitor Suspension Typical Connections on page 185.

Figure 5-18 System with Omega Table and 19" Monitors




(1) The system is provided with one UPS

Figure 5-19 System with Omega Table and Large Display Monitors




- (1) The system is provided with one UPS
- (2) MIS 29414, MIS 29415 and MIS 29416 are connected only with a Third-Party LDM suspension

Figure 5-20 System with Innova IQ OR Table and 19" Monitors



(1) The system is provided with one UPS

Figure 5-21 System with Innova IQ OR Table and Large Display Monitors



(1) The system is provided with one UPS

(2) MIS 29414, MIS 29415 and MIS 29416 are connected only with a Third-Party LDM suspension

5.8 Third-Party Monitor Suspension Typical Connections

Figure 5-22 Typical connections for monoplane with 19" monitors




Figure 5-23 Typical connections for 1 LDM with the 19" backup monitors at the back of the LDM




Figure 5-24 Typical connections for 2 LDM with 19" backup monitors on an additional suspension





Figure 5-25 Connection of the additional in-room 19" monitor



5.9 System Cable Information

5.9.1 Physical Runs

5.9.1.1 Physical Run Synoptic

5.9.1.1.1 System with Omega Table

Figure 5-26 (For Omega Table) Interconnection Length - System with 1 kVA UPS or 8 kVA UPS





Figure 5-27 (For Omega Table) Interconnection Length - System with 20 kVA UPS (Fluoro UPS)



5.9.1.1.2 System with Innova^{IQ} OR Table

Figure 5-28 (For Innova^{IQ} OR Table) Interconnection Length - System with 8 kVA UPS





Figure 5-29 (For Innova^{IQ} OR Table) Interconnection Length - System with 20 KVA UPS (Fluoro UPS)




5.9.1.2 System Core Matrix

NOTICE

The equipment entry point for the PDU, C-FRT & Fluoro UPS is considered to be on the floor, in the middle of the rear face of each cabinet.

For a description of how to use the following cable group schematics, see below:

Figure 5-30 Example of cables schematic



Cable length data is as follows:

- **Cable Length** = the total cable length, connector to connector (example above is 24 meters).
- **X + Y** = used length for connection within system (example above is 3.7 meters).
- **Cable Length - (X + Y)** = available length for conduit run (example above is 20.3 meters).

Figure 5-31 (For System with Omega Table) From Technical Room to Exam Room




Figure 5-32 (For System with InnovalQ OR Table) From Technical Room to Exam Room




Figure 5-33 From Technical Room to Control Room





Figure 5-34 (For System with Omega Table) From Technical Room to Technical Room

Technical Room



Technical Room





Figure 5-35 (For System with InnovalQ OR Table) From Technical Room to Technical Room

Technical Room



Technical Room




Figure 5-36 From Exam Room to Control Room




Figure 5-37 4-6 monitors suspension





Figure 5-38 LDM suspension



5.9.2 MIS (Master Interconnect System)

The system interconnect cables are described in MIS (Master Interconnect System) documents. These documents specify all interconnections between components within the system.

Reference: For specific Vascular system interconnect maps and connection details, refer to the following:

- System MIS Map
- System MIS Charts

General Guidelines

The System introduces a new system interconnect with a star distribution for all cables from the technical area. The cables are shipped on spools to create cable groups. Cable group 1 for Exam room and cable group 2 for control room. The cable group shall be put in place during the same action. The cables are routed in the same duct.

The HV cables could be pulled separately.

5.9.3 Cable Channeling

5.9.3.1 General

High voltage and power cables must be separated from other cables. Use a separate trough in the duct system, or use a separate conduit. Minimize cable length between the MDP and the PDU to reduce voltage regulation problems and wiring costs.

For information about the cables supplied with your system, please refer to [5.9.1 Physical Runs on page 186](#).

5.9.3.2 Conduit

Separate conduits must be used for power and signal wires. These wires must be kept separated from each other.

Using conduit imposes some important considerations when used with this system. Of primary concern, the majority of cables used are pre-terminated. Pre-termination greatly simplifies interconnection but makes cable-pulling difficult because of the added dimensions of the connectors.

Conduit must be large enough to pass the cable and connector through with all other cables already in the conduit. Also, the size of conduit chosen must allow for future growth. There is the possibility of additional cables being added later as the system is developed and options are added.

The use of conduit is recommended for cables running overhead between rooms, especially when a diagonal run provides the shortest cable path.

5.9.3.3 Electrical Ducts

It's important that electrical ducts have separate compartments for power and signal wires. These wires must be kept separated from each other for proper system operation.

Electrical ducts have advantages, when used with a single room or two adjacent rooms. Electrical ducts combine cabling in a neat and functional appearance, with accessibility and room for expansion.

**NOTE**

Mac-lab cables exit behind the table in the Exam Room.

**NOTE**

For **Fast Link** cable and **Private Link** cable (C-FRT Cabinet - AW station), the static operation bending radius must be at least 4 times the outer cable diameter.


It is the responsibility of the site planner to provide the appropriate solution to the table exit (e.g gas box, Clab II, Tram module, connection interface box).

**NOTE**

Specific Recommendations for installation with GE ECG Device such as MacLab, CardioLab or ComboLab:

- TRAM RAC in Exam Room with cable 2016134-106 routed back to Control Room where the other modules & PC are installed
- If no GE Maclab cable 2016134-106 installed between the TRAM (Exam Room) and the Control Room, need to route it so that installation/connection of Physio module can be made in Control Room.

Figure 5-39 Cable Channeling Layout



NOTICE

In some countries, it is forbidden to run electrical cables and water pipes in the same conduit. In this case, two separate conduits are required.

NOTICE

Raceways or cable trays containing electrical conductors shall not contain any pipe, tube or equal for steam, water, air, gas, drainage or any service other than electrical.



NOTE

Only the MEDRAD Mark 7 injector with extension cable requires a separate duct.

**NOTE**

The Physio cable can run in the same conduit as the Bolus cable. In this case, it is required to have a conduit between the table and the physio gases box.

If no conduit available between rear of table and Control Room (no Remote User Interfaces, no MacLab...), need to define proper cable routing or create new conduit as per PIM requirements.

If there is no physio gases box behind the table in the lay out, find a local solution to hide the hole in the floor and the cable exit.

6 Communication Requirements

6.1 Network Requirements

6.1.1 General Information

The system is provided with an internal firewall unit mounted inside the system cabinet and that allows connection to the hospital network for pushing the DICOM images or for service remote access (InSite-RSvP). This firewall is compatible with 10/100/1000 (Gigabit Ethernet) networks.

The C-FRT Cabinet provides an Ethernet RJ45 plug, but the customer is responsible for providing the Ethernet cable between the system and the hospital network.



NOTE


- Connectivity Solutions and pre-installation checklists are available through your local GEHC sales and service representative.
- For InSite-RSvP connections, see information in subsequent sections.

6.1.2 InSite-RSvP Connection Requirements

Service Connectivity for new systems will be based on the InSite-RSvP Platform which allows to configure a direct Internet connection to the RSvP Server (routers/VPN tunnel no more mandatory). Communication with the RSvP server will be outbound only and require using Transport Layer Security (TLS) over TCP port 443. This is commonly known as an HTTPS (HTTP-Secure) connection.

There will be several ways to connect the system to the RSvP Enterprise Server. See below the main options that might not be all available or authorized at your site depending on actual network constraints or local regulations.

Figure 6-1 Connection to the RSvP Enterprise Server - Example



**NOTE**

- The system allows for DNS configuration or proxy server-based connection to the Internet (Option A & B).
- Connection thru a GE Proxy will be possible in the future (Option C & D).
- In the case the customer does not accept the above connection protocol or regulatory reasons prevent using these types of configurations, the local/regional connectivity teams can provide help to connect through SSL/TLS proxy IP over the site-to-site VPN (Option E).

To make the system connectivity operational before the system installation is finished, ensure the connectivity solution is defined as early as possible during the pre-installation process and proper information are exchanged between the customer Network Administrators and GE HealthCare Sales and/or Service representatives.

For all instructions and support, refer to *RSvP Agent Service Manual*. (This document requires specific access right)

IGS RSvP Quick Guides are available through https://iresolve.cloud.ge-healthcare.net/vascular_wiki/index.php?title=Category:Install.

6.1.3 Connection Configuration Parameters

Firstly, the IP addresses for DL and AW PCs have to be requested to the Customer Network Administrators at the time of pre-install to not delay the installation along with:

- A IP address of the hospital Gateway.
- A Subnet mask.
- If additional routers and/or static routes are used by the hospital, those must also be provided.

Regarding the Remote Service connection, the following information will be required from the site depending on their preferred solution:

- the Final system ID for the site (final registration in the CRM/FFA)
- a Domain Name System (DNS) server IP addresses
- or a Proxy server IP or Domain Name and Port
- if a customer wants to only whitelist the specific URLs, the following are the required URLs for RSvP connectivity:
 - <https://insite.gehealthcare.com>
 - <https://as1-insite.gehealthcare.com>
 - <https://as2-insite.gehealthcare.com>
 - <https://gehealthcare-ns.flexnetoperations.com> (for future use only)
 - <https://download.flexnetoperations.com> (for future use only)

**NOTE**

- The System PC (DL) and System Firewall module configurations may differ depending on the final connectivity solution chosen.
- If needed, refer to *InSite® RSvP Agent User OR Technical Reference Manual* for details regarding requesting an RSvP connection setup for the customer.

Refer to Section [6.3 Privacy and Security Configuration on page 206](#) to access the complete list of parameters related to the Privacy and Security configurations that may apply to your site and impact network configurations.

Important Note:

- To configure and verify the RSvP connection, the following accesses will be required:
 - CRM: Contact the regional field support teams to get access to the correct CRM in the region.
 - FFA: request FFA access through MyAccess <https://gehealthcare.saviyntcloud.com/ECMv6/request/requestHome> after completing the required training.
 - RSvP: Through MyAccess <https://gehealthcare.saviyntcloud.com/ECMv6/request/requestHome>, request RSvP access based on region and modality needed after completing the required training.

General Support InSite Connectivity Home Page: https://insiteplus.cloud.gehealthcare.net/#/case_mgmt

- The Training(s) required for access are listed on the linked support central sites.

6.2 DICOM Requirements

The IGS products are DICOM compliant, allowing them to be connected in a network with other DICOM compliant devices for the exchange of images and data.

In some cases, detailed evaluations of the DICOM implementations of devices are needed to ensure interoperability. For this purpose, the DICOM Conformance Statement can be accessed at <https://www.gehealthcare.com/products/interoperability/dicom/xray-mammography-dicom-conformance-statements>, and the IHE Integration Statement can be accessed at <https://www.gehealthcare.com/products/interoperability/ihe/xray-mammography-acquisition-systems-ihe-integration-statements>.

6.3 Privacy and Security Configuration

The Privacy and Security features available with the System require to be configured according to the security policy requested by the hospital.

To ensure the installation is successful and is not delayed because of missing information, it is required to gather all needed information as part of the pre-install process.

The typical parameters are the one listed below. The complete list is provided in Tab "Security Configuration" of the document *IGS System Installation Prerequisites - DOC2024755*. See also Important Notice below.

- **Machine Account**
- **User Authentication**
- **Authorization**
- **Audit Trail**
- **Malware protection**
- **Network Security**
- **Data Transmission and Protection**
- **Other Requirements**

NOTICE

Always refer to the detailed Checklist provided in the document *IGS System Installation Prerequisites - DOC2024755* available from MyWorkshop. Always use the last revision which will contain all mandatory updates.

For details on the new Privacy and Security features available with this machine, refer to the document *Privacy and Security Manual - P/N 5928744-299* available from MyWorkshop.

Support on Privacy and Security can also be found here: <https://www.gehealthcare.com/productsecurity/products>.



www.gehealthcare.com